



Catalog No.	Size
BC-F1520	15 x 20mm
BC-F2030	20 x 30mm



BVitalFORM

Resorbable collagen membrane

Bvital FORM is a strong yet flexible barrier material manufactured from purified porcine peritoneum tissue. It's incredibly easy to handle both dry and wet and can be cut to fit the surgical site exactly, making it ideal for most barrier membrane applications.

Product Features

- Resorbs in 3-4 months
- High mechanical strength
- Highly purified intact porcine peritoneum tissue
- Easy handling and suturable
- Single, uniform layer easy to place correctly
- Permeable to macromolecules and nutrients

Unique Benefits

Make Bvital FORM your membrane of choice and realise the special benefits it brings:

- Greater strength allowing repositioning once placed[^]
- More stable resorption profile[^]
- Lower level of inflammation compared to non-crosslinked collagen[^]

Why **BVital**FORM?


01

Creates a cell barrier to prevent soft tissue ingrowth in the defect site 02

Contains and protects your bone graft material

03

Supports wound healing in guided tissue and bone regeneration procedures



Indications

- Augmentation around implants placed in immediate extraction sockets or delayed extraction sockets
- Localised ridge augmentation procedures
- Alveolar ridge reconstruction for prosthetic treatment
- Filling of bone defects
- Guided bone regeneration in dehiscence defects and guided tissue regeneration procedures in periodontal defects

^ Li, S. T., Yuen, D., Martin, D., & Lee, N. S. (2015). A comparative study of a new porcine collagen membrane to Bio-Gide®. Science, Technology, Innovation, 1-5.

