

INVACARE®MARINER[™] REHAB SHOWER CHAIR



Features

Aluminum frame and stainless steel hardware Fits over most standard and elongated toilet bowls Padded four-position seat Flip-back padded arms that lock down for safety Swing-away front riggings Tool-less adjustable height footrests Folds into a compact unit for storage and transport Removable pan and pan holder Safety strap included 23" treaded urethane tire and 5" caster models







6891w/ 5" Casters

- 1. Padded four-position seat
- 2. Swing-away front riggings with tool-less adjustable height footrests
- 3. Flip-back padded arms that lock down for safety

13"-18"

20.75"x16"x23.5"



6891 - Four 5" locking casters with 18" seat



Invacare Corporation

www.pro.invacare.com

USA

One Invacare Way Elyria, Ohio 44035-4190 (800) 333-6900

Canada

570 Matheson Blvd. E., Unit 8 Mississauga, Ontario L4Z 4G4 Canada (800) 668-5324

Width Width Model # Dimensions Seat Seat-to-Floor Seat-to-Arm Back HxWxD Between Arms Folded Dimensions Height Height Dimensions 6895 w/ 24" Rear Wheels 39"x26.5"x32" 18.5" 18.25"x18.25" 22.75" 15.5"Hx18"W 12" 8' 6891w/ 5" Casters 39"x22.5"x29" 18.5" 8" 18.25"x18.25" 22.75" 8" 15.5"Hx18"W Weight Weight Units Footrest Clearance Over Limited Extensions Toilet HxWxD Capacity Per Carton Warranty 13"-18" 300 lb. 6895 w/ 24" Rear Wheels 20.25"x16"x26.5" 41 lb. 1 3 Years

32 lb.

300 lb.

1

3 Years

©2022

Invacare Corporation. All rights reserved. Trademarks are identified by the symbols ™ and ®. All trademarks are owned by or licensed to Invacare Corporation unless otherwise noted. Specifications are subject to change without notification. Form No. 03-101 Rev. 3/22 220226

