



### FEATURES

- Rear casters help move the trapeze to the desired position
- Chain is adjustable from edge of headboard to the end of boom by sliding to desired position
- Designed to assist larger individuals to change positions while in bed
- Assembly requires a drill, rubber mallet and socket with wrench

## Bariatric Trapeze

Easy assembly and rugged design make the ProBasics Bariatric Trapeze the first choice for heavier users. The 650 lb user weight capacity is the solution to assist larger individuals to change positions in bed. The overhead boom, with an adjustable chain, lets you place the hand grip anywhere from the edge of the headboard to the end of the boom.

### SPECIFICATIONS

|                               |                 |
|-------------------------------|-----------------|
| Actual Product Weight:        | 100 lbs         |
| Maximum User Weight Capacity: | 650 lbs         |
| Outside Legs Width:           | 48"             |
| Assembled Height:             | 73"             |
| Assembled Length:             | 48.5"           |
| Assembled Width:              | 43"             |
| Warranty:                     | 5 years Limited |
| HCPCS Code:                   | E0912*          |

### SHIPPING *(ships in two cartons)*

|  |                   |
|--|-------------------|
| Carton 1 - 70 lbs,<br>45.6" x 3.1" x 29.5" | UPC: 815067072804 |
| Carton 2 - 40 lbs,<br>51.9" x 7.1" x 31.5" | UPC: 815067077021 |

| Item # | Description                 | UOM  |
|--------|-----------------------------|------|
| PBBT   | ProBasics Bariatric Trapeze | each |

\* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.