




**Herbalizer<sup>®</sup>**  
Relief by Design



## SHARE TO WIN

Share **#ReadyToHerbalize** for your  
chance to win a FREE accessory kit



[www.herbalizer.com](http://www.herbalizer.com) [@herbalizeit](https://www.instagram.com/herbalizeit)

## Quick Start Guide

# IMPORTANT

Visit [www.herablizer.com/howitworks](http://www.herablizer.com/howitworks) for  
Herbalizer tutorials and FAQ



## SELECTING TEMPERATURE MODE

To **SWITCH** between **FAHRENHEIT** and **CELSIUS**,  
press and hold both the - and + on the **CONTROL**  
**BAR** for 2 seconds.

## STORING WHIP

Wrap the **WHIP** around the tray **COUNTERCLOCK-**  
**WISE**. Start by placing the long nozzle on its  
**MATCHING ICON**.

## GRINDER CARD

Increase the surface area of herbs by grating them  
before vaporizing. Scrape herbs downward on  
**GRINDER CARD** (just like grating cheese).

## FILLING BOWL

For optimal extraction, fill bowl **HALF WAY** with  
**LOOSELY** ground herbs.

## REMOVING BOWL

Remove **MAGNETIC BOWL** by **TILTING** it forward  
(do not pull up). Once removed, open bowl with  
quarter turn twist.


## CAUTION

This product gets **HOT**, follow usage instructions  
at all times. Avoid touching bottom of bowl and  
bowl seat.


## LIMITED WARRANTY PROTECT YOUR INVESTMENT

The warranty to your 115 volt device is  
**VOID** if used outside the **US** or **Canada**

Visit [www.herbalizer.com/register](http://www.herbalizer.com/register)  
to register your warranty




## AROMATHERAPY MODE

1. Open **ACCESSORY TRAY** and remove the **AROMAPAD** from the front.
2. Apply 2-5 drops of essential oil on the **AROMAPAD** and return it to the front, oil side out.
3. Press  to select **AROMATHERAPY**
4. Set timer using **CONTROL BAR** - + from one minute to one hour
5. Fan speed is automatically set on **HIGH**. For a more **QUIET** experience press the button once.
6. Close the lid and aromatherapy will begin, and when the timer runs out it will automatically shut off.


## USER INTERFACE

### VAPORTHERAPY MODE

1. Press  to select **VAPORTHERAPY**.
2. Remove **MAGNETIC BOWL** by **TILTING IT FORWARD** (do not pull up). Open bowl with quarter turn twist.
3. Fill bowl **HALF WAY** with **LOOSELY** ground herbs

Or, for **CONCENTRATED EXTRACTS** - apply extract to spare **AROMAPAD** and place in **MAGNETIC BOWL** extract-side **UP**.

4. **SET TEMPERATURE** using the **CONTROL BAR** - +



**INTEGRATED SILICONE WHIP** - insert short-end into bowl. Select fan for optional boost.

**SQUEEZE VALVE BALLOON** - attach squeeze valve balloon by inserting on top of bowl, then press fan.

**FREE STYLE** - No whip or balloon, simply press fan to quickly fill environment with vapor.

**Herbalizer**<sup>®</sup>  
Relief by Design