

Good morning Dr. Jones. Can I have a few minutes? I would like to show you a new abdominal binder called QualiBelly Advanced. I think it really can improve our daily work, save time and get your patients more active and comfortable.



That's a lot to ask of an abdominal binder, Mary. What makes you think this one Should be better than what we use now?

OK, the 3 bands have space at each side of the closures before the stitches that connects them. This creates a slit between the bands that allows us to close one band above and one below the drainage tube.



Exactly! That would not only save time but decrease the risk of infection, and save me from making a mess.

Well, even if it's ONE binder, it acts like 3 separate ones. Each can be tightened to a level of compression that's comfortable for the patient.



So compression per band can be set individually, which means the binder can be adjusted to the type of surgery I do?

It's not all, Dr. Jones. It will save even more time for us nurses. You know how we need to be two people to close the traditional abdominal binder on a large patient? Not with this one. One person can easily do it.



I start to see why you like this binder so much, Mary!

Explain that, Mary. Exactly, dr. Jones. Let me show you how we can have different compression levels at the upper, mid, and lower abdomen in this binder. This would increase the comfort of the patient.



Yes, that would certainly be an advantage.
More comfort could help to get patients
to follow our instructions to wear a binder.

I can just open one band, or 2 if needed for good wound access. And since at least one band is kept closed, it's easy for one nurse to close the other 2 when finished, - even on the largest patients! That again saves both time and staff.



Interesting! So, the benefits are: Lower risk of infection, more comfort, saving time - and need of staff. Anything else?

Copyright © Qualiteam Group Ltd. 2020 All rights reserved Yes, I think so. Perhaps even more importantly, we won't need to cut a hole in the binder to get drainage tubes through, - or other devices the patients might have, like a stoma bag.



Yes, I think it will be easier for patients to be more active, and have an earlier bowel movement. Besides adjustable compression, patients can open the top band for comfort when sitting. It's light, aerated and conforms well to the body. Oh, and there is a handy drainage bottle holder that fits on the binder itself without tape or pins!



OK, OK. I get it, Mary. I understand all the benefits for nurses and patients. Let's try it out, and get some hands-on experience with QualiBelly Advanced.