

★ Storage

Store at RT

☞ Contents

- Product Manual
- 10X TBST Buffer, pH 7.4

ALL PRODUCTS SOLD BY GenDEPOT ARE INTENDED FOR RESEARCH USE ONLY UNLESS OTHERWISE INDICATED. THIS PRODUCT IS NOT INTENDED FOR DIAGNOSTIC OR DRUG PURPOSE

★ Introduction

10X TBST solution consisting of 250mM Tris, 1.4 M NaCl 0.5% Tween 20 in high purity RO water, pH 7.4. Solution is 0.2 μ m filtered and dispensed into durable, PET bottles. Dilute to 1X with high purity water for a solution consisting of 25mM Tris, pH 7.4, 140mM NaCl and 0.05% Tween 20. TBST buffer is used as wash buffer and diluent for ELISA, Western blotting and other immunoassays.

★ Formulation

Components	Conc (10X)
TRIS, pH 7.4	250 mM
NaCl	1.4 M
Tween 20	0.5%