




Stress, age, and a busy lifestyle can take their toll—but don't let your skin give you away! ageLOC® LumiSpa® pairs a one-of-a-kind device that uses patented technology with ageLOC® LumiSpa® products. The result? A powerful treatment that helps diminish the signs of aging to keep you looking radiant.

### SYSTEM OVERVIEW

Are you ready for a first-of-its-kind, dual-action skin care system that delivers dramatic skin renewal and deep-cleansing in a single treatment? Invented by Nu Skin and dermatologist-tested, ageLOC LumiSpa is inspired by gene expression studies of skin protein production.

With its soft silicone treatment heads, ageLOC LumiSpa moves at the precise frequency needed to promote skin renewal, revealing smoother, softer skin after just one use.

LumiSpa gently—yet deeply—cleanses away dirt, oil, makeup, pollutants, and toxins through its proprietary oscillating, pore-purifying action. It leaves skin smoother and purer, while reducing the appearance of pores. And the stimulating, massaging LumiSpa experience will leave you looking forward to your next treatment.

Experience brighter, softer, smoother skin and an energizing facial massage. Simply spend two minutes twice a day to enjoy healthier, youthful looking skin in as little as two weeks.

### KEY FEATURES

- **Device**—ageLOC LumiSpa is a waterproof and rechargeable handheld personal skin care device.
- **Treatment Heads**—ageLOC LumiSpa features three treatment heads: GENTLE, NORMAL, and FIRM. Their soft, nonabrasive silicone surfaces contain antimicrobial silver and are more hygienic than brushes commonly used with skin care devices. The sophisticated designs help to maximize product efficacy while making skin more receptive to subsequent products.
- **Treatment Cleansers**—ageLOC LumiSpa features five treatment cleanser options designed for different skin

types: dry, normal to combination, oily, sensitive, and acne prone. These powerful products are formulated using Nu Skin's proprietary ageLOC ingredient blend. Optimally designed for LumiSpa, they are proven to maximize the skin perfecting and cleansing benefits of the device—revealing healthy, bright, youthful looking skin.

### KEY BENEFITS

- Provides instant benefits so skin feels softer and smoother after just one use.
- Delivers seven skin benefits: softness, smoothness, radiance, clarity, purified skin, reduced pore appearance, and visibly improved skin density and volume.
- Benefits intensify over time when consistently used for two minutes twice daily.
- ageLOC LumiSpa uses proprietary Micropulse Oscillation technology, a pore-purifying action that draws out dirt, oil, makeup, pollutants, and toxins and lifts them away without irritating skin while helping to support healthy skin proteins.
- Soft silicone treatment heads are embedded with antimicrobial silver and are more hygienic, less abrasive, and easier to clean than bristle heads commonly found on the market.
- When used in the morning, LumiSpa energizes the skin for a fresher, smoother, more rejuvenated complexion.
- When used at night, LumiSpa deeply cleanses the skin, helping massage and de-stress the skin.
- Removes dead skin cells, promoting dynamic skin renewal and improved skin brightness and smoothness.
- Delivers brighter, healthier, more youthful looking skin in as little as two weeks.
- Gentle enough for twice-daily use.

- Unique sweeping action lifts away loose cells while reducing appearance of pores.
- Goes beyond cleansing to deliver a skin-perfecting treatment through its proprietary oscillating design.
- Doesn't strip skin, but retains skin's natural moisture balance as it cleanses.
- Formulated with Nu Skin's proprietary ageLOC ingredient blend to target the sources of aging and preserve the look of youth.
- ageLOC LumiSpa device is waterproof and can be used in the shower or other wet environment.

### CLINICALLY PROVEN BENEFITS

A 12-week study using ageLOC LumiSpa revealed that individuals who used the product experienced improvements in skin:

- Softness
- Smoothness
- Cleansing
- Brightness
- Radiance
- Clarity
- Luminosity
- Pore Appearance
- Firmness

Participants noticed several of these improvements after a single use.

### USAGE

For optimal results, use ageLOC LumiSpa for two minutes twice daily, as your cleansing step. Morning and night, treat and cleanse with ageLOC LumiSpa treatment cleanser. Then follow with your preferred toner, serum, and moisturizer.

### AGELOC® LUMISPA® TREATMENT CLEANSERS

ageLOC LumiSpa features five treatment cleanser options designed for different skin types: dry, normal to combination, oily, sensitive, and acne-prone. These powerful products contain Nu Skin's proprietary ageLOC ingredient blend to target the sources of aging and are formulated to maintain healthy, bright, youthful looking skin.

ageLOC treatment cleansers are specifically formulated to enable the precise cushioning, cleansing, and interaction with the skin, optimizing the effectiveness of LumiSpa's proprietary treatment head action.

### AGELOC® LUMISPA® TREATMENT CLEANSERS KEY INGREDIENTS

- **ageLOC LumiSpa Treatment Cleanser—Dry** contains Sodium Lactate, Sodium PCA, Glycine, Fructose, Niacinamide, and Inositol to help protect skin's natural moisture, Squalane to help maintain skin's natural moisture barrier, and Salicyloyl Phytosphingosine to soothe the skin.
- **ageLOC LumiSpa Treatment Cleanser—Normal/Combo** contains Milk Thistle (*Silybum Marianum* Seed) Oil, *Porphyra Umbilicalis* Extract, and *Nannochloropsis Oculata* Extract to help maintain skin's natural moisture barrier, as well as *Rhododendron (Ferrugineum)* Extract to provide antioxidant protection.
- **ageLOC LumiSpa Treatment Cleanser—Oily** contains Japanese Knotweed (*Polygonum Cuspidatum* Root) Extract, Myristyl Alcohol, and PCA to help reduce oil production and soothe the skin, as well as Pomegranate (*Punica Granatum*) Juice Extract and Carnosine to provide antioxidant protection.
- **ageLOC LumiSpa Treatment Cleanser—Sensitive** is our most gentle formula. It is a pink gel cream formulated with a mild surfactant system, which results in a low foam. It contains *Avena Sativa (Oat) Kernel* Extract, *Pinus Tabulaeformis* Bark Extract (which contributes to the pink color), and Bisabolol to help calm and soothe the delicate skin. Allantoin also helps to condition the skin.
- **ageLOC LumiSpa Treatment Cleanser—Acne** is a mild formula that contains 0.5% salicylic acid to help clear pores and reduce most breakouts. It is low-foaming to prevent additional skin irritation. It also contains Carnosine to provide antioxidant protection and *Epilobium Angustifolium* Flower/Leaf/Stem Extract that can help calm and support acne-prone skin.

### COMPLEMENTARY PRODUCTS

- **ageLOC Galvanic Spa®/Nu Skin® Facial Spa**—ageLOC Galvanic Spa uses direct/galvanic currents to rejuvenate your complexion and increase skin vibrancy. Nu Skin Facial Spa uses microcurrents to stimulate and tone the skin.
- **ageLOC Me®**—creates an entirely new approach to customized anti-aging skin care. Brilliant yet simple to use, ageLOC Me delivers five powerful anti-aging products, custom selected from almost 2,000 possible regimen combinations. Delivering one streamlined, simplified regimen to you.
- **Tru Face® Line Corrector**—targets signs of aging with peptides that help soften the appearance of lines around your mouth, eyes, and forehead.

# ageLOC® LumiSpa®

- ageLOC® Tru Face® Essence Ultra—a powerful, lightweight serum that targets aging at its source to visibly firm and contour skin while protecting against damaging free radicals.
- ageLOC® LumiSpa® Accent and IdealEyes®—a LumiSpa attachment and eye treatment duo that harnesses the technology of ageLOC LumiSpa to improve the appearance of the eye area.

## FREQUENTLY ASKED QUESTIONS

### How does ageLOC LumiSpa work?

With its soft silicone treatment head, ageLOC LumiSpa's proprietary Micropulse Oscillation technology enables the device to move at the precise frequency needed to promote skin renewal, revealing smoother, softer skin after just one use. At the same time, LumiSpa gently—yet deeply—cleanses away dirt, oil, makeup, pollutants, and toxins through its proprietary oscillating, pore-purifying action, leaving skin smoother and purer while reducing the appearance of pores.

### What are the main differentiators of ageLOC LumiSpa?

There are three key differentiators of ageLOC LumiSpa:

1. Unlike cleansing devices currently on the market, which feature abrasive brush or bristle heads, ageLOC LumiSpa features a gentle, more hygienic silicone treatment head that won't damage the skin or lift skin cells.
2. ageLOC LumiSpa features proprietary Micropulse Oscillation technology to stimulate the skin to more effectively treat and cleanse the skin and deliver superior skin renewal benefits.
3. ageLOC LumiSpa is a dual-action, timesaving device that delivers a superior skin-perfecting treatment and deep cleansing benefits simultaneously.

### What is the rationale behind the name?

"Lumi" stands for luminous and emphasizes the initial benefits that using the device and products will have on your skin. "Spa" reinforces the fact that the device, treatment heads, and treatment cleansers will deliver spa-quality results every time you use them.

### How do I know which ageLOC LumiSpa treatment head is right for me?

Select your ageLOC LumiSpa treatment head based on your skin type and treatment needs. Each treatment head removes impurities and purifies skin, delivering a gentle experience ideal for a variety of skin types. The Gentle treatment head is our gentlest choice and is great for those with sensitive or acne-prone skin. The Normal treatment head is the standard choice and is designed

for most skin types. And the Firm treatment head is for those who want the most aggressive experience but is not recommended for those with sensitive or acne-prone skin. You can use different heads for different experiences each day.

### How long will my treatment head last?

We recommend replacing your treatment head every three months to ensure your skin will receive optimal benefits on a consistent basis.

### How does ageLOC LumiSpa work with the treatment cleansers?

The treatment cleansers were developed specifically for use with ageLOC LumiSpa. They're designed to work synergistically with the oscillating motion of the heads to optimize distribution of anti-aging ingredients onto the skin while cleansing and exfoliating to promote healthy, youthful looking skin.

### How long should the treatment cleanser last?

One treatment cleanser should last 30 days.

### Are there any patents on ageLOC LumiSpa?

Yes, there are multiple patents issued and pending for ageLOC LumiSpa and the proprietary Micropulse Oscillation technology.

### Do I need to use the ageLOC LumiSpa treatment cleansers with ageLOC LumiSpa, or can I use other treatment formulas or cleansers with it?

ageLOC LumiSpa treatment cleansers are specifically formulated for use with the ageLOC LumiSpa device. They enable the precise cushioning, cleansing, and interaction with the skin, optimizing the effectiveness of LumiSpa's proprietary treatment head action. They are the only formulations tested and approved for use with ageLOC LumiSpa to deliver the full system benefits and promote healthy, youthful looking skin.

### How is ageLOC LumiSpa different from ageLOC Galvanic Spa® or Nu Skin® Facial Spa?

ageLOC LumiSpa provides twice daily treatment and cleansing and should be used as part of your regular skin care routine, whereas ageLOC Galvanic Spa and Nu Skin Facial Spa are specialized treatments that can be used three times per week. ageLOC LumiSpa and ageLOC Galvanic Spa/Nu Skin Facial Spa deliver complementary anti-aging benefits: ageLOC LumiSpa treats

# ageLOC® LumiSpa®

and cleanses skin for a smoother, softer, more radiant appearance while ageLOC Galvanic Spa conditions, revives, and hydrates the skin, and Nu Skin Facial Spa uses microcurrents to tone and contour the skin.

## Does ageLOC LumiSpa treat acne?

ageLOC LumiSpa was not designed to treat acne; however, ageLOC LumiSpa Treatment Cleanser—Acne is a mild formula designed for acne-prone skin. Because it contains 0.5% salicylic acid, the cleanser helps clear pores and reduce most breakouts. You should use the Gentle head with ageLOC LumiSpa Treatment Cleanser—Acne. The ageLOC LumiSpa system was developed to promote healthy, youthful looking skin and to deliver general skin renewal and purification benefits. Those treating severe acne should continue with their regular skin care regimen and consult a dermatologist before using ageLOC LumiSpa.

## Do I need to wash my makeup off before using LumiSpa?

ageLOC LumiSpa conveniently and effectively cleanses away makeup, so you do not need to wash makeup off before using the device. However, to remove heavy makeup, you might want to consider using another cleanser in addition to your LumiSpa. LumiSpa should not be used on the eyes, so you should use a separate eye makeup remover to remove makeup on and around the eyes.

## How do I know if I'm using the right motion and right amount of pressure?

ageLOC LumiSpa has a built-in system that will guide you in using the right motion and right amount of pressure. If you use a scrubbing motion or press too hard, the device will pause and vibrate until you apply normal speed and pressure. Slowly move the device across your skin in broad strokes. The device will pause briefly every 30 seconds to prompt you to move from one area of the face to the next.

## Can I use ageLOC LumiSpa on my body?

ageLOC LumiSpa was designed to target facial skin and deliver anti-aging and deep cleansing benefits in this area. We do not recommend it be used elsewhere on the body.

## Can I use ageLOC LumiSpa in the shower?

Yes, ageLOC LumiSpa is waterproof, and you can use it in the shower or other wet environments.

## How long does ageLOC LumiSpa last after a full charge?

ageLOC LumiSpa was developed to last at least a week between charges with regular use.

## INGREDIENTS

### AGELOC LUMISPA TREATMENT CLEANSER—DRY

Water (Aqua), Dicaprylyl Carbonate, Diisostearyl Malate, Tribehenin PEG-20 Esters, Glycerin, Glycereth-26, Squalane, Cetearyl Alcohol, Dimethicone, Caprylic/Capric/Myristic/Stearic Triglyceride, Sodium Astrocaryum Murumuru, Oryza Sativa (Rice) Extract, Schizandra Chinensis Fruit Extract, Xanthophyll, Narcissus Tazetta Bulb Extract, Acrylates/C10-30 Alkyl Acrylate Crosspolymer, Niacinamide, Urea, Methylsilanol Hydroxyproline Aspartate, PEG-40 Hydrogenated Castor Oil, Sodium Lactate, Sodium PCA, Vegetable Oil, Xanthan Gum, Batyl Alcohol, Butylene Glycol, Pentylene Glycol, Phytosteryl/Octyldodecyl Lauroyl Glutamate, Maltodextrin, Brassicamidopropyl Dimethylamine, Hydroxyphenyl Propamidobenzoic Acid, Salicyloyl Phytosphingosine, Fructose, Glycine, Inositol, Lactic Acid, Lecithin, Phenoxyethanol, Hydroxyacetophenone, Disodium EDTA, Sodium Hydroxide, Fragrance (Parfum).

### AGELOC LUMISPA TREATMENT CLEANSER—NORMAL/COMBO

Water (Aqua), Coco-Betaine, Sodium Cocoyl Glycinate, Glycerin, Acrylates Copolymer, Diisopropyl Sebacate, Glycol Distearate, PEG-8, Rhododendron Ferrugineum Extract, Nannochloropsis Oculata Extract, Schizandra Chinensis Fruit Extract, Porphyra Umbilicalis Extract, Silybum Marianum Seed Oil, Narcissus Tazetta Bulb Extract, Oryza Sativa (Rice) Extract, Acrylates/C10-30 Alkyl Acrylate Crosspolymer, Polyquaternium-39, Phenoxyethanol, Chlorphenesin, Sodium Benzoate, Phytic Acid, Sodium Hydroxide, Sodium Phytate, Fragrance (Parfum), Synthetic Fluorphlogopite, Titanium Dioxide (CI 77891), Alumina, Tin Oxide.

### AGELOC LUMISPA TREATMENT CLEANSER—OILY

Water (Aqua), Cocamidopropyl Betaine, Acrylates Copolymer, Ammonium Laureth Sulfate, Glycerin, PEG-80 Sorbitan Laurate, Di-PPG-2 Myreth-10 Adipate, Coco-Glucoside, Punica Granatum Juice Extract, Narcissus Tazetta Bulb Extract, Polygonum Cuspidatum Root Extract, Schizandra Chinensis Fruit Extract, Oryza Sativa (Rice) Extract, Glyceryl Oleate, Polyquaternium-10, Myristyl Alcohol, PCA, Carnosine, 1,2-Hexanediol, Fragrance (Parfum), Citric Acid, Aminomethyl Propanol,

# ageLOC® LumiSpa®

Dehydroacetic Acid, Chlorphenesin, Phenoxyethanol, Sodium Benzoate, Propanediol, Caprylyl Glycol, Disodium EDTA, Mica (CI 77019), Titanium Dioxide (CI 77891).

#### **AGELOC LUMISPA TREATMENT CLEANSER—SENSITIVE**

Water (Aqua), Cetearyl Alcohol, Glycerin, PEG-8 Dimethicone, Cocamidopropyl Betaine, Caprylyl/Capryl Glucoside, Capryl/Capramidopropyl Betaine, Carbomer, Narcissus Tazetta Bulb Extract, Pinus Tabulaeformis Bark Extract, Schizandra Chinensis Fruit Extract, Avena Sativa (Oat) Kernel Extract, Oryza Sativa (Rice) Extract, Allantoin, Zingiber Officinale (Ginger) Root Extract, Bisabolol, Carnosine, Xanthan Gum, Hydroxyacetophenone, Phenoxyethanol, Disodium EDTA, Sodium Hydroxide, Sodium Chloride, Citric Acid.

#### **AGELOC LUMISPA TREATMENT CLEANSER—ACNE**

Water (Aqua), Cetearyl Alcohol, Glycerin, PEG-8 Dimethicone, Cocamidopropyl Betaine, Caprylyl/Capryl Glucoside, Capryl/Capramidopropyl Betaine, Carbomer, Salicylic Acid, Narcissus Tazetta Bulb Extract, Epilobium Angustifolium Flower/Leaf/Stem Extract, Schizandra Chinensis Fruit Extract, Xanthan Gum, Carnosine, Oryza Sativa (Rice) Extract, Disodium EDTA, Phenoxyethanol, Hydroxyacetophenone, Sodium Hydroxide, Sodium Chloride.

