

WCS™ Closed-Toe

Wound Care Shoe System



The WCS™ Closed-Toe Wound Care Shoe System is the product of choice when dealing with Wagner Grade 3 or greater ulcerations. The four multi-density insoles can be modified, effectively decreasing pressure from the plantar aspect of the foot.



DARCO



Replacement insoles available. Sold separately.

WCS™ Closed-Toe Features and Benefits:

- › **Four Multi-density Insoles** allow for the ultimate in customization and precisely targeted off-loading.
- › **Ultra-soft Plastazote® Lining** provides a virtually friction-free contact surface that won't exacerbate existing wounds.
- › **Sold in Pairs** to eliminate discrepancies in shoe height and the potential associated knee and hip pain.

Insoles Cross Section

- Plastazote®
Thickness 1/4" Durometer 26 →
- Medium Density EVA
Thickness 1/4" Durometer 42 →
- Firm Density EVA
Thickness 1/4" Durometer 50 →
- PORON® Cushioning
Thickness 1/8" Durometer 29 →



WCS™ Closed-Toe:

- › Wagner Grade 3 or greater ulcerations
- › Bony prominences
- › Long-term treatment of wounds and ulcerations
- › Maximum off-loading while still allowing for patient ambulation

| Item Number | Size | Color | Women's Size | Men's Size |
|-------------|--------------|-------|--------------|------------|
| WCC0B | X-Small | Black | 4 - 5.5 | - |
| WCC1B | Small | Black | 6 - 8.5 | 6 - 7.5 |
| WCC7B | Small/Medium | Black | 9 - 9.5 | 8 - 8.5 |
| WCC2B | Medium | Black | 10 - 10.5 | 9 - 9.5 |
| WCC3B | Large | Black | 11 - 12.5 | 10 - 11.5 |
| WCC4B | X-Large | Black | 13 - 14.5 | 12 - 13.5 |
| WCC5B | 2X-Large | Black | 15 - 17 | 14 - 16 |

Available through:



Headquarters:

DARCO International | 810 Memorial Blvd. | Huntington, WV 25701 | USA

Toll Free voice 800-999-8866 | Fax 877-201-7082 | Voice 304-522-4883 | Fax 304-522-0037

info@darcointernational.com | www.darcointernational.com