

Validated 3D Printers For pro3dure Resins Carried By Gro3X



3D Printer – material settings **Dental**



		3D Printer								
available models		ASIGA MAX	ASIGA MAX UV	ASIGA PRO2	Dazz 3D L120	MiiCraft Ultra	MiiCraft Prime	Structo Velox	UNIZ SLASH PLUS	phrozen Shuffle/Sonic
Generative Resin	GR-10 guide	✓	✓	✓	✗	✓	○	✓	✗	✗
	GR-11 tray	✓	✓	✓	✗	✓	○	○	✗	✗
	GR-12 cast	✓	✓	✓	✗	✓	○	○	✗	✗
	GR-12 LCD cast	✗	✗	✗	○	✗	✗	✗	✓	○
	GR-13 model	✓	✓	✓	✗	○	✓	○	✗	✗
	GR-13 LCD model	✗	✗	✗	✓	✗	✗	✗	✓	○
	GR-15 gingiva	✓	✓	✓	✗	○	○	○	✗	✗
	GR-16 Xray	✓	✓	✓	✗	✓	○	○	✗	✗
	GR-17 temporary	✓	✓	✓	✓	✓	○	○	○	○
	GR-17.1 temporary It	✓	✓	✓	✓	✓	○	✓	○	○
	GR-18.1 IB	✓	✓	✓	✗	✓	○	✓	✗	✗
	GR-19 OA	✓	✓	✓	✗	○	○	✓	✗	✗

pro3dure medical works continuously on the development and updating of parameters for selected open 3D printing systems.

- ✓ available
- in progress
- ✗ not available

Materials are our DNA!

pro3dure medical GmbH
 Am Burgberg 13 · 58642 Iserlohn, Germany
 Phone: +49 (0)2374 920050-0
 Fax +49 (0)2374 920050-50

info@pro3dure.com · www.pro3dure.com