

On The Go 27 pcs/box Art.No: 110 **1-K** Ready to use

FreeBact® Water

Drink 50L 27 pcs/box Art.No: 106

Tank 50-250L 27 pcs/box Art.No: 141



Tank 250-500L 27 pcs/box Art.No: 142


Tank 1-2 m³ 27 pcs/box Art.No: 143

Tank 2,5-5 m³ 12 pcs/box Art.No: 145

Tank 5-10 m³ 12 pcs/box Art.No: 147



XINIX FreeBact[®] EMERGENCY PRODUCTS Military & Civil Safe drinking water and surfaces



XINIX

EEBACT

XINIX

- 99,8% water and 0,2% chlorine dioxide (ClO2).
- 65-95% cheaper than bottled water.
- Minimal plastic waste.
- Does not contain chlorine or alcohol.
- Eliminates the biofilm in the tank.
- 5 times more efficient than chlorine/3 times more than alcohol.
- 1-K products activated within minutes
- (compare to chlorine tablets: 30 mins waiting time).
- Cleans the water and the tank at the same time.
- Non flammable i.e. possible to transport by air.
- You drink no chemicals, only barely measurable amounts of salt.
- By Swedish Food Agency approved drinking water additive.
- Produced in Sweden by a Swedish company, founded 1992.
- No aftertaste or smell.
- Safe & gentle for both people, animals & nature.
- More environmentally friendly than chlorine.
- Both customers, partners and lab-tests agree; Xinix products are superior to chlorine based products.
- Scientifically proven.
- Used today by leading companies, authorities and aid organizations worldwide.
- Easy to use (see other side).





Xinix Emergency Box

56 bottles - pocket size → á 50-250 liters

→ This box will generate 14.000 liters of safe water!





1-K Ready to use



FACTS

Chlorine dioxide (ClO₂) is a GAS by nature – but Xinix has managed to stabilize and bind it into LIQUID FORM.

Chlorine dioxide is not a new disinfectant; it was discovered 1811 by Sir Humphry Davy. It is well-renowned as one of the worlds most efficient disinfectants. It has since then primarily been used for treatment of drinking water, but on a bigger scale only in industries. With Xinix FreeBact® now available in your pocket!

After treating your water with Xinix FreeBact[®], the only residue is barely measurable amounts of salt, NaCl. You can't even taste it. Compare this to other water disinfectants, where you actually might drink chlorine. With Xinix FreeBact[®] all the good minerals will also remain in the water.

FreeBact[®] SURFACE On The Go



2 years **Opened bottle**: 30 days Still active when yellow

Perfect 50 ml spray to keep in your bag or jacket. Effective & easy-to-use disinfection for all types of surfaces. Spray your hands, handles, tables, fruit - and rest assured that you are not bringing unwanted bacteria or viruses with you.



OTHER WATER / with visible particles FreeBact[®] WATER TAP WATER / clear water Drink 02 03 02 01 100 ml 1 min 50 Liters 100 ml 2 min With this little bottle (50 ml) in the bag you get safe water. To be dropped into your glass or mug. Drink after 1-2 min. Shelf life: years Add 100 ml of water, pour through a coffee filter, a tissue or similar Opened bottle Add 100 ml of water Wait 1 minute before drinking Wait 2 minutes before drinking 5 drops into vessel l drop into a 30 days drinking glas Still active when uellow. Work in an open space or ventilated area. FreeBact[®] WATER Clean tank & water: 6 Tank (different sizes) \cap 1 box of Xinix FreeBact Box 1: step 1-4 + 7-8 Unscrew the two childproof 50 - 10.000 Liters caps on the bottles. В A Clean tank, water & pipes: 2 boxes of Xinix FreeBac Box 1: step 1-6 Clean up the tank and get rid of 02 \rightarrow Read more: Add all of solution B (small neck) Box 2: 1-4 + 7-8 unwanted biofilm – at the same time В <u>xinix.se</u> to solution A (large neck). as uou disinfect the drinking water. 05 No need to refill with new water. Warning: A strong odor will Pour the solution into the tank occur. Never inhale the gas Drink after 40 min. 06 before or after mixing. Shelf life: Flush all pipes. Mixed bottles: Fill tank with water again 24 hou 03 Screw the cap back on the A-bottle and shake it gently. Ð 07 Pour the solution into Place on a steady surface. 30 min the tank. Wait 30 min before use 04 10 min Open the cap half a turn. 08 Wait 10 minutes before using Э Your tank and water (+ piping system) the mixed solution. are now disinfected.



