


## SRx Mini-Pearls™ Compounding System

Excellent for making a high volume of small dose troches or suppositories quickly. The Mini-Pearls System has 228 small spherical cavities with a diameter of 6.3mm and a volume of 0.2ml ea. (228 x 0.2ml = 45.6 ml plate total). Average cavity weight is 0.2g of Techna 20 Troche Base. Durable, flexible, non-stick white silicone top mold, that is oven, microwave, freezer, and dishwasher safe. Temp resistant range -60 to 230°C (-76 to 446°F). Lightweight yet rigid high-density black polymer set plate will not scratch countertops and has no sharp edges. Includes four stainless steel alignment pegs.



**\*\*Included in Kit:**

DESCRIPTION	PART #	COUNT	U/M
228 Cavity White Silicone Mini-Pearls Top Mold	E0591-00	1	EA
228 Cavity Black Set Plate, 4 x SS Posts Inserted	E0589-00	1	EA
Roller	E0594-00	1	EA
Scraper Spatula	E0569-00	1	EA
60ml Catheter Tip Syringe	E0568-00	1	EA
Mini Hot Air Heating Tool	E0570-00	1	EA
Tote, Laboratory	E0592-00	1	EA
Tote Cover, Laboratory	E0593-00	1	EA
Instruction Manual- Use and Cleaning		1	EA
Techna 20 SF Troche Base (Sample)	B0070-08	1	EA

\*\*See our website or call us at 888-512-1209 for current kit and component pricing.