# SPECIALIZED PRESCRIPTION COMPOUNDING SYSTEM SRx Mini-Pearls<sup>TM</sup> System

# For Compounding Small Dose Troches & Suppositories

- High volume troches or suppositories quickly
- 228 X 0.1 ml spherical cavities
- Optimal shape for patient comfort
- Quicker, easier, smaller

Watch our SRx Mini-Pearls<sup>™</sup> System tutorial video at: https://www.youtube.com/watch?v=xF7yUREVyGU

### SRx Mini-Pearls System Kit (E0590-00) includes:

- Black 228 Cavity Set Plate
- White Silicone 228 Cavity Top Mold
- Mini Hot Air Heating Tool
- Mini-Pearls Roller
- Scraper Spatula
- 10ml Oral Tip Syringe
- Tote, Laboratory Multi-Purpose
- Techna 20 SF™ Troche Base Sample
- Step-by-step Use and Cleaning Guide

## SRx Mini-Pearls<sup>™</sup> Features

The SRx Mini-Pearls System has 228 small spherical cavities with a diameter of 6.3mm and a volume of 0.13ml ea. (228 x 0.13ml = 29.64 ml plate total). Average Mini-Pearl weight is 0.2g of Techna 20 SF™ Troche Base.

The non-stick white silicone top mold is durable & flexible, as well as oven, microwave, freezer, and dishwasher safe: temp resistant range -60 to 230°C (-76 to 446°F). The white color enables easy visual cleanliness inspection.

Lightweight yet rigid black high-density polymer set plate will not scratch countertops and has no sharp edges: temp resistant range -50°C to 90°C (-58°F to 194°F) and includes four stainless steel alignment pegs.

To watch system instructional video go to: https://www.youtube.com/watch?v=xF7yUREVyGU

SRx Mini-Pearls™ System Kit Item#E0590-00 \$595.00

Additional 228 Cavity Silicone Top Mold Item#E0591-00 ......\$80.00

#### **Troches Benefits**

- Unique small spherical shaped dose for patient comfort
- Increased patient compliance
- Fits up to 80mg of API per cavity
- Quick dissolution
- Increased buccal absorption of active drug
- Faster onset of action
- Reduction of dental caries

#### **Suppository Benefits**

- Unique small spherical shaped dose for patient comfort
- Easy administration vs creams/gels
- Mini-Pearls fit perfectly in most vaginal/ rectal applicators
- Fits up to 80mg of API per cavity
- Quick dissolution
- Faster absorption of active drug
- Faster onset of action

## Order Today!

# Call 1.888.512.1209 www.specializedrx.com

email: sales@specializedrx.com



SpecializedRx Products, LLC | 7237 University Ave NE | Fridley, MN 55432 | USA



For bulk orders, please call us.