

# ALCOHOL FREE enduro HAND SANITISER

99.99%  
kill rate against  
Enveloped  
Viruses  
EN 14476  
TESTED & PROVEN EFFECTIVE



Quick drying, long lasting, water-based alcohol free Enduro Hand Sanitiser is designed for everyday use - anytime, anywhere! Dermatologically tested, Enduro Hand Sanitiser dries within seconds, leaving hands clean, sanitised and naturally moisturised.

## Overview

- **Effective against enveloped viruses**
  - tested according to International Standard EN 14476
- **Kills up to 99.9999% of bacteria & fungi**
  - starts working in seconds
- **Long lasting protection**
  - proven efficacy up to 4 hours after application
- **Naturally moisturises**
  - leaves skin clean and refreshed
- **Dermatologically tested**
  - suitable for repeated everyday use
- **Alcohol free**
  - with no odour or fragrance and no sticky residue

## How to use



Pump one or two applications into palm of hand.



Rub gently into front and back of hands until dry.



Don't forget finger nails, skin creases and between fingers.

## When & where to use



On the move



Hospitals, clinics & care homes



Eating out



At work



Kitchens & food prep areas



Travelling



Commuting



Before and after visiting the toilet



Offshore

## Test results

Tested in accordance with the following European standards against:

<b>Viruses*</b>	EN 14476
<b>Bacteria</b>	EN 1276 EN 14561
<b>Fungi</b>	EN 1275 EN 1650
<b>Spores</b>	EN 13704

### Other results

Dermatologically tested  
Passes 96 hour skin patch

30 second kill rates  
Tested against bacteria

2 hour residuality  
Tested against bacteria, fungi & spores

4 hour residuality  
Tested against bacteria

Full details on our test data is available to professional users on request: to make an enquiry, please contact us.

\* Enveloped viruses: limited virucidal spectrum activity

USE BIOCIDES SAFELY. ALWAYS READ THE LABEL AND PRODUCT INFORMATION BEFORE USE.



Bio Technics Ltd.  
Linton Business Park, Gourdon, Aberdeenshire, Scotland, UK DD10 0NH  
Tel: +44 (0) 1561 361515 Email: info@biotechnics.co.uk Web: www.biotechnics.co.uk



FM 696382 EMS 696323

# Pack sizes



## Dispenser bottles

Description	Case size	Order code
50ml foam pump dispenser	Case of 12	HS-50F
100ml foam pump dispenser	Case of 20	HS-100F
150ml foam pump dispenser	Case of 20	HS-150F
200ml foam pump dispenser	Case of 16	HS-200F
250ml foam pump dispenser	Case of 16	HS-250F
500ml foam pump dispenser	Case of 12	HS-500F
750ml foaming trigger sprays	Case of 12	HS-750TS



## Dispenser refills

Description	Case size	Order code
1 Litre bottles	Case of 6	HS-1L-RF6
1 Litre bottles	Case of 12	HS-1L-RF
5 Litre container	Case of 2	HS-5L-RF2
5 Litre container	Case of 4	HS-5L-RF



## Wall-mounted dispenser & refills

Description	Capacity	Case size	Order code
Automatic dispenser - pouch fill	1 Litre	Sold individually	AUTW-1000P
1 Litre pouch refill for AUTW-1000P		Case of 12	AUTW-1000P-HS
Institutional manual dispenser - pouch fill	1 Litre	Sold individually	MANW-INST
1 Litre pouch refill for MANW-INST		Case of 12	MANW-INST-HS
Manual dispenser - pouch fill	0.8 Litre	Sold individually	MANW-800P
0.8 Litre pouch refill for MANW-800P		Case of 12	MANW-800P-HS
Manual dispenser - reservoir fill	1 Litre	Sold individually	MANW-1000R



## Bulk containers

Description	Case size	Order code
20 Litre container	Sold individually	HS-20
200 Litre drum	Sold individually	HS-200
1000 Litre IBC	Sold individually	HS-1000

## Also available as a concentrate (1:4 dilution)

Description	Case size	Order code
200 Litre drum	Sold individually	HS-C-200
1000 Litre IBC	Sold individually	HS-C-1000

Manufactured by



**Bio Technics Ltd.**

Linton Business Park, Gourdon, Aberdeenshire, Scotland,

Tel: +44 (0) 1561 361515 info@bio-technics.com



**SYNERGY** MEDICAL SUPPLY  
BRINGING QUALITY AND SERVICE, WE CARE!

Tel: (852) 2892 0033 Fax: (852) 2892 1433  
E-mail: admin@synemed.info