



alkoCLEAN

gauze pad moistened with 70% isopropyl alcohol

Revision I
Date 2016.06.08

Page 1/1

MANUFACTURER : Yinjing Medical Supplies Co. Ltd.
DISTRIBUTOR: ZARYS International Group sp. z o.o. sp.k

CONTACT : ul. Pod Borem 18, 41-808 Zabrze, Poland, tel. +48 32 376 07 00-02

INDICATION:

- skin disinfection before injection, vaccination or blood sampling

FEATURES:

- made of non-woven
- soaked with isopropyl alcohol (70%), which is a very good disinfectant and has bactericidal and fungicidal effect
- two pouches connected in one-piece packaging can be separated thanks to perforation



SIZE AND REFERENCE:

REF	DIMENSIONS	NO. OF PLIES	PACKAGING UNIT (PAPER-ALUMINIUM POUCH)	INTERMEDIATE PACKAGING	BULK PACKAGING
L-02	30 mm x 65 mm (unfolded) 30 mm x 32,5 mm (folded)	2	2 x 1 pc.	100 x 1 pc. (50 x 2 pcs)	100 boxes
L-03	56 mm x 65 mm (unfolded) 28 mm x 32,5 mm (folded)	4	2 x 1 pc.	100 x 1 pc. (50 x 2 pcs)	100 boxes



Zarys International Group sp. z o.o. sp.k.,
ul. Pod Borem 18, 41-808 Zabrze, Poland
tel. +48 32 271 69 91 fax. +48 32 274 72 84
www.zarys.com.pl