

MANUFACTURER: ZARYS International Group sp. z o.o. sp.k.


ADDRESS: ul. Pod Borem 18, 41-808 Zabrze, Poland, phone +48 32 376 07 00-02, fax. +48 32 370 38 08

INDICATION:

- Device designed for managing the flow of infusion fluids during administration
- Allows to switch the infusion flow path
- For blood pressure monitoring
- For direct drug / fluid supply

FEATURES:

- Three-way stopcock with an extension tube of various lengths (from 7cm to 100cm)
- Transparent design
- Rotating Luer Lock
- Clear optical identifier of open/close position (marked flow directions)
- Made of polycarbonate
- Drain made of PVC
- All inlets protected with stoppers
- Resistant to lipids
- Latex-free
- DEHP-free
- Single-use
- Non-pyrogenic, non-toxic
- EO sterilized
- Shelf life: 5 years
- Individual packaging: paper-foil



REFERENCE

CATALOGUE NUMBER	LENGTH OF THE TUBE	INDIVIDUAL PACKAGING	INTERMEDIATE PACKAGING	BULK PACKAGING
KTP-7	7 cm	paper-foil	50 pcs	20 x 50 pcs
KTP-10	10 cm	paper-foil	50 pcs	20 x 50 pcs
KTP-25	25 cm	paper-foil	50 pcs	20 x 50 pcs
KTP-50	50 cm	paper-foil	50 pcs	20 x 50 pcs
KTP-75	75 cm	paper-foil	50 pcs	20 x 50 pcs
KTP-100	100 cm	paper-foil	50 pcs	20 x 50 pcs