



Advanced antimicrobial defense for your surfaces

SIQURA™ Hospital Grade Surface Disinfectant & Protectant destroys bacteria and locks onto surfaces to provide 24hrs of active defense on high touch surfaces and up to 30 days on low touch areas. It kills 99.999% of bacteria and is effective against viruses* including COVID-19.

Most disinfectant products evaporate within a short period of time, leaving surfaces vulnerable to recontamination. **SIQURA™ 75HG** is a three-in-one solution. It **cleans, disinfects** and then remains bonded to create a **self-sanitising** surface that continues to defend against bacterial attack.

24 up to 24 hrs protection on high touch surfaces

30 days up to 30 days protection on low touch surfaces

kills 99.999% bacteria and viruses*


bonds to most surfaces

easy-to-use rapid action

How do you use SIQURA™ 75HG?



1 Ensure the surface is clean and dry. Spray and wipe **SIQURA™ 75HG** over the entire surface making sure it is completely covered.



2 Allow to air-dry. Do not wipe surface as this will remove the coating before it's had a chance to bond. Once dry, the antibacterial active is locked onto the surface.



3 Surface is now secured and will defend against bacterial contamination for up to 24hrs on high touch areas and 30 days on low touch areas, unless heavy abrasion occurs.

How friendly is SIQURA™ 75HG?

Rather than poisoning microbes, **SIQURA™ 75HG** uses a small electrostatic charge to kill microorganisms, making it suitable for use in healthcare, childcare, education, commercial and domestic environments. It can be used on a wide range of surfaces and high touch areas including counters, door handles, bathrooms, floors, furniture and seating. Rapidly degrading, SIQURA™ 75HG breaks down to form simple organic compounds — reducing the potential for negative environmental impacts.

- ✓ Solvent-free
- ✓ No chlorines
- ✓ No bleaches
- ✓ No carcinogenic compounds
- ✓ Non-oxidative
- ✓ Non-flammable
- ✓ Bio-friendly

How durable is it?

SIQURA™ 75HG uses Quatlock technology to bond the active antimicrobial to most surfaces for long-lasting defense. For more information on the science behind SIQURA, please visit siqura.com.au. **SIQURA™ 75HG** is listed on the ARTG (Australian Register of Therapeutic Goods) as a Hospital Grade Surface Disinfectant that protects surfaces for up to 24 hours in high touch areas and up to 30 days in low touch areas.

What is it effective against?

***SIQURA™ 75HG** has been tested according to the following methods EN1276, EN13697, EN13624, ASTM 1053, EN 14476, AOAC 991.47, AOAC 991.48, AOAC 991.49 and ASTM E2149 and is effective against a broad range of bacteria and viruses such as Human Coronavirus E229 (including COVID-19), Murine Hepatitis Virus (MHV), Staphylococcus aureus, Enterococcus hirer, Escherichia coli, Pseudomonas aeruginosa and Salmonella choleraesuis.

How is it available?

CODE	DESCRIPTION	PACK
0-1528	75HG 500ml spray head bottle	12
0-1634	75HG 5 litre drum	2
0-1573	75HG 20 litre drum	1

Ideal for use within:

- Healthcare:** Hospitals & Medical Clinics
- Hospitality:** Restaurants & Cafes
- Education:** Childcare, Schools, Universities
- Public Transport:** Buses, Trains, Trams
- Retail Spaces:** Counters, Touchpoints
- Lifestyle Environments:** Gyms & Fitness
- Homes:** Furniture, Floors, Fixtures
- Vehicles:** Cars, Vans, Trucks

