Astroskin **Product Specifications**



Physiological Data Signals















3-lead ECG

256 Hz, 6.5uV Resolution

- Heart rate: 30-220 BPM, 1 Hz
- QRS event detection: 4 ms resolution RR intervals: 4 ms resolution

Quality assessment channels:

- Disconnection detection
- 50-60 Hz noise detection
- Saturated ECG signal detection
- Movement artifact detection RR intervals reliability detection

Dual Channel Breathing Sensors

(RIP) 128 Hz

- Breathing Rate: 3-60 BPM, 1 Hz
- Minute Ventilation: 80-65000 mL, 1 Hz
- Tidal Volume (last inspiration): 80-10000 mL, 1 Hz, 20 mL resolution
- Inspiration and Expiration Events: 8 ms resolution

Quality assessment channels:

- Disconnection detection Noise detection
- Baseline change detection

Pulse Oximetry

- Oxygen saturation (SpO2%): 1Hz
- Photoplethysmography (PPG): 64 Hz
- Heart rate: 1 Hz

Systolic Blood Pressure (BP)

- 1 Hz
- Pulse Transit Time computation

3-axis Accelerometer

64Hz, +/-16g range, 0.004g resolution.

- Actigraphy: 1 Hz
- Step counting: reported at each step detection: 30-240 SPM
- Cadence: 1 Hz (0, 30-240 SPM, calculated on the 8 last steps)

Skin Temperature

1Hz, 0.1 Celsius resolution





Astroskin Shirt - For Women



Astroskin **Product Specifications**



OF CONSECUTIVE RECORDING

Physical Specifications

Physical Overview

- Dimensions: 85 x 55 x 25 mm

Device weight: 89 gGarment weight: < 100 g

Battery

2x AA batteriesAutonomy: 48 h

Embedded Local Data Storage

- Capacity: 4 GB of data

- Autonomy: ~ 560 h of recording

MADE IN CANADA

This Garment is Made with Italian fabric

Composition:

73% Polyamide Micro 27% Elastane



Excellent UV Protection



Antibacterial Treatment



Exceptional Breathability



Piling and Streching Resistant



Machine washable



Quick Drying

Astroskin Research Platform

- iOS Apps iPhone & iPad Real-time Data Visualization
- Data Synchronization
- Dashboard
- API



INSERT


48 HOURS

2 AA BATTERIES



Astroskin Shirt Sizing

With built-in all-textile sensors, the Astroskin Shirt is light and comfortable.



User Interface Button

Battery Compartment Door

Women Size _s o	2XS	xs	S	M	L	XL	XXL
Thorax	Inches 26 - 27.5	28 - 29.5	30 - 31.5	32 - 33.5	34 - 35.5	36 - 37.5	38 - 40
	Centimeters 65 - 70	70 - 75	75 - 80	80 - 85	85 - 90	85 - 90	90 - 95
Abdomen	Inches 26 - 27.5	28 - 29.5	30 - 31.5	32 - 33.5	34 - 35.5	36 - 38	38.5 - 40.5
	Centimeters 65 - 70	70 - 75	75 - 80	80 - 85	85 - 90	90 - 95	95 - 100

