Caremed Alrick.

Mobilise CH5 chair

An innovative early mobilisation chair*, made in the UK



An innovative, early mobilisation chair that lays 100% flat and has an innovative legrest that encourages sit to standing, this chair is manufactured in the UK, using UK materials which greatly reduces leadtimes.

Manual Handling

Lay flat option reduces the need for hoists



Retractable legrest

Fully retracts beneath the seat so easier sit to standing



Easy to move

Allows patients to be moved outside or to other wards if necessary



Versatile

Simple and easy to use handset with icons for ease of use



Seating options

Depending on the needs of the patient, different seating options are available



The Mobilise chair is made in the UK and therefore has shorter leadtimes



Features

- Multi-purpose early mobilisation chair
- ✓ Lays 100% flat
- Sit to standing option
- 160kg safe working load
- Retractable legrest fully retracts beneath the seat so easier sit to standing
- Seat tilt with Trendelenburg
- Twin wheel 125mm castors with total lock or directional lock
- Adjustable seat width 510-560mm
- Easily remove the seat and back cushions for cleaning
- Thick cushion base with different seating options depending on patient
- Tough, soft touch Dartex®

 Endurance material, with minimal stitching
- Strong, padded removable armrests
- Lay flat dimensions; 1830mm length by 560mm wide

Contact us for a demo or purchase:

updates@caremed.healthcare 01926 889677