


Scan the QR code to watch instruction video



Cushion Size Guide

Step 1

Print the measuring tape on an A4 or 8 1/2" x 11" paper.
Please confirm that your print is scaled correctly by placing a credit card in the indicated outline. It should fit perfectly.



Step 2

Cut out the measuring tape.
Place the measuring tape on the inside of your closed watch bracelet (not on your wrist).



Step 3

Determine your cushion size. If you are between cushion sizes, please reach out to info@charles-simon.com and a member of our team will assist you in finding the right fit for you.

- Small** 14.5 - 16 cm / 5.75 - 6.25 in
- Medium** 16 - 17.25 cm / 6.25 - 6.75 in
- Large** 17.25 - 18.5 cm / 6.75 - 7.25 in

