

300 Technology Drive  
Christiansburg, VA 24073 USA  
inorganicventures.com

P: 800-669-6799/540-585-3030  
F: 540-585-3012  
info@inorganicventures.com

## **ACCREDITATION / REGISTRATION**

Inorganic Ventures is an ISO 9001 registered manufacturer (QSR Certificate Number QSR-1034).

## **DISCRETE ANALYZER REAGENT** Ascorbic acid solution

Catalog No: DA-PO4-ASCORBIC  
Lot No: T2-DA728143


Starting Material: Ascorbic acid  
Starting Material Lot No: IV-2590

**NOMINAL CONCENTRATION:** 15 g/L Ascorbic acid

**This solution is a reagent and is not intended to be used as a certified reference material.**

## **PREPARATION**

This reagent solution is ready to use without further preparation for analyzing o-Phosphate in water utilizing discrete analyzer methods based on method chemistry described in EPA Method 365.1. Figure 1 demonstrates a o-Phosphate calibration curve generated on a SEAL AQ400 discrete analyzer using this lot of reagent in combination with the additional IV method reagent(s): DA-PO4-COLOR



**Figure 1:** Correlation Coefficient = 0.9997 / Polynomial Order = 1

Concentration	Result	Acceptance Limits
0.02	0.0095	0.0056 – 0.0112
0.05	0.0209	0.0122 – 0.0245
0.10	0.0400	0.0236 – 0.0473
0.25	0.0966	0.0608 – 0.1213
0.50	0.1956	0.1276 – 0.2126
0.75	0.3008	0.2017 – 0.3361
1.00	0.4042	0.2635 – 0.4391

**STORAGE AND HANDLING**

Store refrigerated at 4° C while in sealed TCT bag. After opening the sealed TCT bag keep cap tightly sealed when not in use. Continue to store bottle at 4° C when not in use. Do not return removed aliquots to container. Allow for reagent to reach room temperature and shake well before use.

Certification Date: 01/04/2023

Prepared By:

*Uyen Truong*

Approved By:

*Michael Z Booth*

Expiration Date: 06/04/2023

Certifying Officer:

*Paul R. Gaines*