

EMPOWERING WOMEN'S WELLNESS WITH INMODE'S APPLICATOR VERSATILITY (SURGICAL OR NON-SURGICAL)



AVIVA

Included on Surgical Platform

Aviva is a state-of-the art, **minimally invasive technology**, designed for electro-coagulation of the Fibro Septal Network (FSN), sub-dermal fat, dermal collagen remodeling, and hemostasis. Designed with extensive safety features, Aviva reduces the risk of thermal injury during the remodeling of the soft tissue matrix e.g. non-surgical labiaplasty.

Treatment areas include:

- Labia Majora
- Mons Pubis
- Introitus
- Labia Minora
- Perineal Body



Specifications

| | |
|-----------------------|---|
| Cannula Length: | 6cm |
| Cannula Diameter: | 0.9mm |
| Internal Temperature: | Operator adjustable between 50°C – 65°C |
| External Temperature: | Operator adjustable up to 38°C |



TONE


Can be included on Surgical or Non Surgical Platform

Tone is a non-invasive muscle strengthening technology that works by emitting low frequency pulses delivering **electrical muscle stimulation (EMS)** to induce involuntary muscle contractions. Patients rest comfortably while Tone targets specific muscle groups to **restore abdominal strength** without surgery or downtime.



Specifications

| | |
|----------------------------|-----------------------|
| Intensity (output voltage) | Up to 50 energy level |
| Pulse Width | 20-400 [usec] |



inmode.com.au/pages/empower-rf-womens-wellness



[InMode Australia](#)
[InMode_Australia](#)
[InModeAustralia](#)
[盈美特 InMode Australia](#)
[458006879](#)

INFO@INMODE.COM.AU | INMODE.COM.AU | +61 (2) 7209 5678



EMPOWER RF
by INMODE



Imagine
the
Difference

EMPOWERING WOMEN'S WELLNESS





THE NEW STANDARD EMPOWERING WOMEN'S WELLNESS

EmpowerRF is the next generation women's wellness technology platform, delivering life changing treatments. EmpowerRF leverages three complimentary energies:

- submucosal fractional radiofrequency (Fractional RF micro-needling)
- intravaginal electrical muscle stimulation (EMS)
- bipolar radiofrequency

This intelligent programmable device comes equipped standard with four synergistic technology hand pieces and two expandable applicators:

- Morpheus8V
- Morpheus8
- VTone
- Aviva
- FormaV
- Tone

The unique combination of multiple modalities allows for customised, layered procedures to address a broad range of age groups and feminine conditions including:

- Urinary incontinence:
 - Stress (SUI)
 - Urge (also known as over active bladder)
 - Mixed
- Sexual dysfunction
- Pelvic floor weakness
- Genitourinary syndrome of menopause (GSM) or vulvovaginal atrophy (VVA)
- Changes in the aesthetic appearance of the female anatomy
- Dryness, itching, reduced elasticity, and burning during menopause

EmpowerRF helps to restore women's self-confidence and quality of life without surgery or prolonged downtime.


"EmpowerRF is a game changer for women's health. It's the platform with all technologies that complement each other to improve muscle tone, rejuvenate the intravaginal area and treat urinary incontinence".

DR ELIZABETH GOLEZ
Cosmetic Physician & Obstetrician
(Caring for women's health for over 25 years)



MORPHEUS8 V

Morpheus8V reinforces and tightens the vaginal stroma and basal lamina, connective tissue trophism, and helps reduce SUI through Morpheus8's clinically proven triple action of electrocoagulation, collagenous tissue contraction, and sub-necrotic bulk heating. With a 24 pin matrix of thin gold-coated microneedles and a treatment depth of 3mm, Morpheus8V delivers targeted and precise energy delivery, achieving surgical-like outcomes without general anaesthesia or downtime associated with excisional interventions.



Specifications

| | |
|------------------|---|
| Configuration: | 24 gold coated micropins |
| Depth Profile: | 1mm - 3mm |
| Pin Profile: | 300um diameter with insulated 0.5 mm conductive tip |
| Treatment Depth: | Operator adjustable |



V TONE

VTone is a gentle, nonsurgical technology designed to provide intravaginal electrical muscle stimulation (EMS) and neuromuscular re-education to rehabilitate weak pelvic floor muscles. The female anatomy evolves over time due to multiple factors: pregnancy, childbirth, hormonal transformations, menopause, and aging. Inevitably, the female pelvic floor is negatively impacted leading to pelvic floor weakness, stress urinary incontinence, and a reduction in the quality of life. VTone's two-sided tip inserts quickly and comfortably into the vaginal canal and delivers the maximum contractions to treat vulvovaginal atrophy. It is a highly efficacious and painless non-invasive procedure to improve SUI symptoms.




Specifications

| | |
|-----------------------|---|
| Applicator: | Designed to comfortably treat the female anatomy |
| Technology: | Electrical Muscle Stimulation (EMS) - 72,000 pulses/ 30 min treatment |
| Treatment Parameters: | Operator adjustable |



MORPHEUS8

Morpheus8 is a fractional radiofrequency microneedling technology providing a triple action of sub-necrotic bulk heating, coagulation of fat and connective tissue, and ablation. RF energy is delivered through bipolar arrays of pin electrodes, and results in localised heating and coagulation of the tissue that is in direct contact with the pins. The result is improvement in skin surface irregularities, skin tone, laxity, fine lines, wrinkles and indented scars. Fractional coagulation promotes remodeling while untreated tissue between the pins promotes faster healing. Morpheus8 can positively impact health and quality of life.



Tip Specifications:

| | RESURFACING | PRIME | MORPHEUS8 | MORPHEUS8 body |
|------------------|--------------|---------------------|---------------------|---------------------|
| Configuration: | 24 Micropins | 12 Micropins | 24 Micropins | 40 Micropins |
| Depth Profile: | 0.5mm | 1mm - 4mm | 1mm - 4mm | 2 - 7mm |
| Treatment Depth: | Fixed | Operator adjustable | Operator adjustable | Operator adjustable |

FORMA V

FormaV helps to reduce/improve dryness and sexual dysfunction by vaginal tightening, due to total metabolic thermal stimulation, increased blood flow, proliferation of the mucosa, hydration, and physiological trans-tissue activation of the vagina walls.

FormaV is a comfortable treatment that provides uniform volumetric deep heating for tissue remodeling with real-time temperature control.



Specifications

| | |
|----------------------|--------------------------------|
| RF Output Power: | Up to 40W |
| Temperature Cut-off: | Operator adjustable up to 43°C |

