

Catheter Slides

PART #	DIMENSIONS (D inches)	SPECIFYING IN LOGICELL
CAT16HD	For 16 shelf/cabinet	Allow space for catheters to hang
CAT24HD	For 24 shelf/cabinet	Allow space for catheters to hang



CAT16HD
Catheter Slide

CAPACITY OR MAXIMUM WEIGHT LOAD

CAT16HD: There are 6 hooks, each hook can hold up to 15 packages.
CAT24HD: There are 9 hooks, each hook can hold up to 15 packages.

Label Holders

DESCRIPTION

These plastic label holders attach to the front of trays and or baskets to provide a location for labeling the contents of the unit for improved organization and access.

USE

The LH3 is attached by a strap which is wrapped around a horizontal wire of the basket and therefore can be placed in various locations across the front of the basket. The LH132PS and LH134PS are mounted by sliding on the full width of the top edge of the basket or tray. The clear labeling area of the LH132/134 is folded over on itself to provide a surface for placing the label as well as a cover to protect the label.

Label Holders

PART #	DIMENSIONS (H x L inches)	DESCRIPTION
LH3	1" h x 3"	Specify color: white, yellow, orange, pink, purple, green, red, or blue*
LH132PS	2" h x 13"	For trays and baskets
LH134PS	4" h x 13"	For trays and baskets

* For baskets only.



LH134PS
Label Holder



LH3RED
Label Holder



LH132PS
Label Holder