



**1624SHC**Shelf with Cath Rack Channel



244DIV Divider



1624SHC & 244DIV
Shelf with Divider Assembly

# **Shelves and Dividers**

## **DESCRIPTION**

These shelves are used to store larger or bulky items in LogiCell units. The shelf has a 2.5"h lip on the rear to prevent items from sliding back, and the front drops down 1"h to provide a labeling surface, but shelf can be inserted in the reversed direction to retain items on shelf in slanted position.

#### USE

Shelf can be slanted. Items ending in "C" have an added channel across the rear and front that accepts catheter slides (2416SHC up to 3, 1624SHC up to 5, 2436SHC up to 7). Slots every 1" for attaching dividers.

# **Shelves and Dividers**

PART#	<b>DIMENSIONS</b> (D x W x H inches)	SPECIFYING IN LOGICELL	DIVIDER
1624SHC or 1624SHE	16 x 24 x 2.5	For a 24"w column. With divider, 2 cell slots used	164DIV (14" long) maximum of 21
2416SHC or 2416SHE	23.5 x 15 7/8 x 2.5	For a 16"w column. With divider, 2 cell slots used	244DIV (21 1/2" long) maximum of 13
2436SHC or 2436SHE	23.5 x 32.75 x 2.5	For a 33"w column. With divider, 2 cell slots used	244DIV (21 1/2" long) maximum of 31

#### **MATERIALS**

Shelf and dividers are made of 18g steel, epoxy coated, white finish.

## CAPACITY OR MAXIMUM WEIGHT LOAD

Items ending in "C": 70# Items ending in "E": 35#

# **Catheter Slides**

## **DESCRIPTION**

These LogiCell slides are used to store catheters contained in pouches. The rack pulls-out of the cabinet for ease of access. The front plate has a label holder (3" x 5") for improved organization.

# **USE**

The slide can be attached to the ceiling of a LogiCell cart / cabinet, or below the xxxxSHC metal shelf.