Creating Bespoke Solutions Together





Amazing Skin [Activation Kit]



Bespoke Solutions for:

- anti-aging
- radiance
- skin brightening

treatments of the:

- face
- neck
- décolleté

FREQUENCY

Recommended course of 6 to 10 sessions, as follows:

1 session/week for 1 month

1 session every 2 weeks for 1 month

1 session after 6 weeks

1 session every 3 to 6 months





HYALURONIC ACID (3.5%)

Volumizing Solution

Box of 5 vials - 5ml

- Volume
- Wrinkle Eraser
- Water Binder



HYALURONIDASE (0.5%)

Depolymerizing Enzyme

Treating reticulated Hyaluronic Acid granulomas

Box of 5 vials - 5ml

How to use:

1 vial + 4ml Saline Solution

0.1 - 0.2 cc per point









DMAE (5%)

Firming Solution

Box of 20 ampoules - 5ml

- Lifting
- Flaccidity (tummy, arms, legs)
- Sagging skin (body & face



SILICIUM (1%)

Draining & Regenerating

Box of 20 ampoules - 5 ml

- Anti ageing
- Skin Rejuvenation
- Wound healing
- Alopecia
- Lipolysis
- Cellulite
- Stretch Marks











GLUTATHIONE (5.5%)

Lightening Solution

Box of 5 vials

- Photo-ageing
- Melasma



ASCORBIC ACID (20%)

Brightening Solution

Box of 5 vials - 5 ml

- Photo-ageing
- Freckles
- Dull Skin
- Smoker's skin









GLYCOLIC ACID (4.5%)

Rejuvenating Solution

Box of 20 ampoules - 2 ml

- Dry skin, Sun-damaged skin
- Fine Lines and Wrinkles
- Loss of Elasticity
- Large Pores, Bumps
- Age Spots / Dark Spots
- Actinic Keratosis



RETINOL 300 (15%)

Revitalizing / Resurfacing Solution

Box of 20 ampoules - 2 ml

- Sun damaged skin
- Smoker's skin
- Fine Lines
- Acne
- Wound healing
- Stretch Marks









Protocols - Face



Foundation for Regeneration Skin Imperfection/ Recovered Wound		Severe Imperfection	
		Product needed	Per Session
Product needed	Per Session	(SPN6)	4.5.2.1
SPN-5	3ml	SPN-5	1.5 -3ml
		Collagen Booster	2ml

Comp	lementary	v Combin	ations
Comp	icincincal)	Combin	utions

Rejuvenation & Brightening

Product needed	Per Session
Skin Renew	1.5ml
Glutathione	1.5ml

Firming & Radiance

Per Session
1.5ml
1.5ml

Acne Scars & Open Pores

Product needed	Per Session
Glycolic Acid	1ml
Hyaluronic Acid 3.5	1ml
Retinol 300	1ml

Tailored Skin Anti-Ageing Mix

Product needed	Per Session
Skin Renew/Skin Lift	1ml
Glycolic Acid	1ml
DMAE	1ml

Protocols - Eye & Neck



Foundation for Regeneration

Skin Imperfection/ Recovered Wound

Product needed	Per Session
SPN-5	3ml

Severe	Imperfectio	n

Product needed (SPN6)	Per Session
SPN-5	1.5 -3ml
Collagen Booster	2ml

Complementary Combinations

Eye Area (mild microneedling)

Product needed	Per Session
Skin Lift	0.4ml
DMAE	0.3ml
HA 3.5	0.3ml

Neck Treatment

Product needed	Per Session
Skin Renew	0.5ml
DMAE	1ml
Hyaluronic Acid 3.5	1ml
Ascorbic Acid	0.5ml
Retinol (optional)	1ml