

Salmon Range Protocol

SPN-6

Supporting Skin Regeneration

For rejuvenation and

skin imperfections



SPN-6 | Package content:

- 1 vial of 3ml SPN-5 Solution
- 1 vial of 2ml Collagen Booster


> Supporting skin regeneration at a glance

INDICATIONS & MAIN INGREDIENTS:

SPN-6 (SPN-5 + Collagen Booster) | Long-lasting skin regeneration for deep wrinkles and severe skin imperfections. Indicated for the treatment of the face, neck and body.

Collagen Booster Main Ingredients: Polycaprolactone, Hyaluronic Acid.

> Recommended protocol | SPN-6





Supporting Skin Regeneration

For rejuvenation and

skin imperfections



SPN-5 | Package content:5 vials of 3ml SPN 5 Solution

> Supporting skin regeneration at a glance

INDICATIONS & MAIN INGREDIENTS:

SPN-5 | Supporting skin regeneration, rejuvenation, wrinkles and light skin imperfections. Indicated for the treatment of the face, neck and décolleté.

SPN-5 Main Ingredients: PDRN, Hyaluronic Acid, Glutathione, Acetyl Hexapeptide-8, Palmitoyl Pentapeptide-4.

> Recommended protocol | SPN-5

Light skin imperfections Use 3ml SPN-5 per treatment area Skin texture and brightening Use 5ml SPN-5 per treatment area

> Recommended treatment frequency: 5-7 sessions

