

LIVEN C^o

Ring sizing guide

We always recommend having your ring size measured in-person at a jewelry store. In the meantime, please reference the below methods to obtain your approximate ring size.



Print this page, ensuring that it is set to print at 100% without scaling. To make sure the scale is correct, measure the 2 lines below. They should be 1 inch, and 50mm.

1 inch

50mm


Option one: Measure an existing ring

If you have a ring that is a good fit, you can use the guide below to find out its approximate size. Simply place your ring over the circle that fits without any gray showing inside - remember, you are measuring the *internal diameter* of the ring. For additional guidance, please email info@livenco.com, or call us at (213) 293-6088.



Option two: Measure your finger

Again, ensure your printer is set to 100% with no scaling, and cut out the below. Make a slit where the line is, and thread the end through the slit around your finger. Make sure you can slide the sizer over your knuckle. We recommend measuring your finger 3-4 times, as your size can change slightly due to the time of day, temperature etc.



Please note, we are not able to guarantee ring sizing via a paper method. Please confirm your size via a jewelry professional if possible.