

# Thermo Scientific Arctic Express Cryogenic Shippers with Datalogger

Maintain the chain and your peace of mind.



Maintaining the cold chain is paramount to your samples, including when they are in transit. Thermo Scientific™ Arctic Express™ cryogenic shippers are designed to ensure the maintenance of sample integrity when movement of your samples is required.

Arctic Express models are now available with a built-in 21CFR, part 11 compliant datalogger, providing you with additional peace of mind that your important samples have been maintained during transport.

#### Shipping containers

- Innovative material absorbs liquid nitrogen to prevent spillage during shipment
- 14- and 21-day liquid nitrogen holding times allow specimens to be shipped safely

- Durable construction and sturdy base allow dry shipper to withstand rough handling
- Lockable lid prevents unauthorized entry
- Easily transported via most common carriers
- Lightweight aluminum design and convenient pail-style handle

#### Datalogger

- 21CFR11 compliant
- High- and low-programmable alarm settings
- Temperature precision +/- 0.5°C
- Digital display
- Easy to use: automatically generate PDF report, summarizing shipping temperatures



**Not available outside North America**

**Thermo**  
SCIENTIFIC

## Thermo Scientific Arctic Express Cryogenic Shippers with Datalogger (Not available outside North America)

	Arctic Express 5	Arctic Express 10	Arctic Express 10 IATA	Arctic Express 20
<b>Model No.</b>	<b>CY50915D</b>	<b>CY50905D</b>	<b>CY50905DIATA</b>	<b>CY50910D</b>
LN2 capacity (liters)	1.5	4.3	4.3	10
Canisters per unit	1	1	1	N/A
Canes per canister	3	8	2	161
2mL vials per cane	3	6	6	6
2mL vial capacity	9	48	12	966* or 500**
Total straw capacity	15	120	20	1610
Static holding time, <-150°C (days)	8	21	21	14
Neck diameter in. (cm)	1.4 (3.5)	2 (5.0) 8.5 (21.6)	2 (5.0) 8.5 (21.6)	8.5 (21.6)
Exterior dimensions Diameter x height: in. (cm)	7.3 x 13.5 (18.5 x 34.3)	8.7 x 19.4 (22.1 x 49.3)	8.7 x 19.4 (22.1 x 49.3)	15 x 23 (38.1 x 58.4)
Primary padlock	Included	Included	Included	Included
Secondary padlock	Included	Included	Included	Included
Shipping shell	Included	Included	Included	Included
Datalogger (21CFR11)	Included	Included	Included	Included
Datalogger replacement cap	N/A	Option (see below)	Option (see below)	Option (see below)
IATA approved Kit	N/A	N/A	Included***	N/A

\* 966 bulk, six vials per cane

\*\* 500 vials using five 2-inch boxes with 100-cell dividers

\*\*\* Included with initial purchase of Arctic Express 10 IATA (CY50905DIATA). IATA Kit is also available as accessory and should be used with Arctic Express 10 products (CY50905, CY50905D and CY50905DIATA) only.

### Optional Accessories

Model No.	Description
AEDL10	Cork/cap assembly with datalogger for Arctic Express 10
AEDL20	Cork/cap assembly with datalogger for Arctic Express 20
AEDLACC	Additional datalogger for Arctic Express units
GT261910	IATA approved packout kit* for Arctic Express 10. NOTE: ONLY for use with models CY50905, CY50905D or CY50905DIATA

\* Kit contains IATA tested Primary Receptacle, Secondary Packaging and Rigid Outer Packaging

[www.thermoscientific.com/cold](http://www.thermoscientific.com/cold)

© 2015 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries. Specifications, terms and pricing are subject to change. Not all products are available in all countries. Please consult your local sales representative for details.

**Australia:** +61 39757 4300

**Austria:** +43 1 801 40 0

**Belgium:** +32 53 73 42 41

**China:** +86 21 6865 4588 or  
+86 10 8419 3588

**France:** +33 2 2803 2180

**Germany national toll-free:** 0800 1 536 376

**Germany international:** +49 6184 90 6000

**India toll free:** 1800 22 8374

**India:** +91 22 6716 2200

**Italy:** +39 02 95059552

**Japan:** +81 45 453 9220

**Netherlands:** +31 76 579 55 55

**New Zealand:** +64 9 980 6700

**Nordic/Baltic/CIS Countries:**

+358 9 329 10200

**Russia:** +7 812 703 42 15

**Spain/Portugal:** +34 93 223 09 18

**Switzerland:** +41 44 454 12 22

**UK/Ireland:** +44 870 609 9203

**USA/Canada:** +1 866 984 3766 (866-9-THERMO)

**Other Asian Countries:** +852 2885 4613

**Countries not listed:** +49 6184 90 6000

**Thermo**  
SCIENTIFIC

Part of Thermo Fisher Scientific