

You can have it all in a MaxiComfort® lift and recline chair. Brought to you exclusively by the world's leading motion experts at Golden Technologies.



SIT Sit 

The **SIT** button moves the recliner to the proper default seated position automatically. MaxiComfort offers an exclusive true 90° back angle found easily by using the arrows at the bottom of the hand control. MaxiComfort also features a back recline without raising the footrest.



ZG Zero-gravity 

The **ZG** position was created to eliminate pressure on the body. A body in the zero-gravity position will decompress after a long day by gliding into pressure free comfort and experience a natural relaxation with minimum muscle tension.




TV TV Watching 

The **TV** Watching position was designed to offer the most comfortable, neutral posture for upright relaxation. A body in the TV Watching position is recharged with ergonomic upper body support and comfort.

SLEEP

With both the back and the feet in the horizontal "layflat position," enjoy the recliner for a restful sleep.



Trendelenburg 

The therapeutic Trendelenburg position was created to improve circulation in the body. Once in the ZG position, the recliner can further elevate the legs and feet above the heart to increase blood flow and release tension from the lower extremities. A body in the Trendelenburg position receives an elongation of the lumbar area for a rejuvenating stretch.

