

# AirLife® ETCO<sub>2</sub> sampling cannulas

AirLife® ETCO<sub>2</sub> sampling cannulas offer a soft, clear design to optimize patient comfort while allowing unobstructed patient observation. Available in pediatric or adult configurations to meet hospital ETCO<sub>2</sub> monitoring needs, these cannulas:

- Incorporate a permanent barrier separating the nasal prongs.
- Allow simultaneous O<sub>2</sub> inhalation and end-tidal gas sampling.
- Help deliver O<sub>2</sub> to one naris while sampling CO<sub>2</sub> from the other, for more accurate ETCO<sub>2</sub> readings.
- Offer crush-resistant tubing to prevent occlusion and kinking, adding to patient safety.
- Are made without DEHP.

Cat. no.	Description	Qty. (cs.)
With male luer connectors		
2802M-10	Pediatric with 7' O <sub>2</sub> tubing, 7' CO <sub>2</sub> tubing	10
2811M-10	Adult with 7' O <sub>2</sub> tubing, 2" CO <sub>2</sub> tubing	10
2812M-10	Adult with 7' O <sub>2</sub> tubing, 7' CO <sub>2</sub> tubing	10
2812M-10U	Adult with 7' U/Connect-It® O <sub>2</sub> tubing, 7' CO <sub>2</sub> tubing	10
2812M-25	Adult with 7' O <sub>2</sub> tubing, 7' CO <sub>2</sub> tubing	25
With female luer connectors		
2802F-10	Pediatric with 7' O <sub>2</sub> tubing, 7' CO <sub>2</sub> tubing	10
2811F-10	Adult with 7' O <sub>2</sub> tubing, 2" CO <sub>2</sub> tubing	10
2812F-10	Adult with 7' O <sub>2</sub> tubing, 7' CO <sub>2</sub> tubing	10
2812F-25	Adult with 7' O <sub>2</sub> tubing, 7' CO <sub>2</sub> tubing	25
ETCO <sub>2</sub> sampling lines		
28M2M	ETCO <sub>2</sub> sampling line 10', male-to-male luer connectors	10
28M2F	ETCO <sub>2</sub> sampling line 10', male-to-female luer connectors	10



To learn more, contact your CareFusion sales representative at **800.323.9088** or visit [carefusion.com/AirLife](http://carefusion.com/AirLife)

**⚠ WARNING**—U.S. Federal Law restricts this device to sale by or on the order of a physician.

CareFusion  
Vernon Hills, IL

[carefusion.com](http://carefusion.com)

