

- > 5-10CM distance measurement
- Avoid cross infection

Imported sensors

- > measurement deviation ≤ ±0.2
- > 1 second measurement

Three color backlight

- Convenient night measuremen, more direct indication for body temperature
- Automatic shutdown after 30s of no operation

- Low fever: yellow backlight
 - + beep to remind
- > High fever: red backlight + beep to remind

32 sets of memory

- Grasp the temperature change
- Support clinical thermal diagnosis

Body temperature and temperature dual mode

- Body temperature, milk temperature, water temperature, room temperature, baby care in all directions
- Safe and healthy, dual temperature unit change between degrees Celsius and Fahrenheit