

# Cocoon N95 Disposable Respirator Masks

# Specifications

- Particle Filtration Efficiency >95% for 0.3 Micron
- Synthetic Blood Penetration at 160mmHg (level 3)
- Compliant to N95 NIOSH standard (42 CFR 82 subpart K), tested by Intertek USA and Nelson Lab USA
- Compliant to AS 1716, tested by Flinders University, South Australia
- Made in Australia
- ARTG 340518

#### **Order Codes:**

MSK 001 Cocoon N95 Mask -standard size

MSK 002 Cocoon N95 Mask with Headband

MSK 003 Cocoon Disposable N95 Mask with Headband and foam nose wire

Pack size: 20 units per box

#### **MSK 001**



**MSK 002** 



**MSK 003** 





Head office: 25 Slevin Street, North Geelong Victoria 3215, Australia

Factory: 103 Mornington Street, North Geelong Victoria 3215, Australia Care Essentials Pty Ltd ABN: 89 148 060 486

Telephone: +61 3 52771455 Fax: +61 3 52771411



#### **Order Codes:**

**MSK 001S** Cocoon Disposable N95 Mask – small size

**MSK001SF** Cocoon Disposable N95 Mask - small size with foam nose wire

Pack size: 40 units per box

#### **MSK 001S**



**MSK 001SF** 



## **Mask Accessories**

#### **Ear Savers**

- Made in Australia using FDA approved polypropylene using injection molding technology.
- Can be used with both surgical and N95 masks.
- Helps achieve good fit and reduces strain on ears.

### **Order Code:**

ES 001 Ear Saver







Head office: 25 Slevin Street, North Geelong Victoria 3215, Australia

Factory: 103 Mornington Street, North Geelong Victoria 3215, Australia Care Essentials Pty Ltd ABN: 89 148 060 486

Telephone: +61 3 52771455 Fax: +61 3 52771411