

# Male Catheterization/Enema Simulator

MW52



# Female Catheterization/Enema Simulator II





## practice Achieve better training results with our soft and

High quality simulators for urethral catheterization, perineal care and enema

realistic catheterization simulator that can also be used on SP or full body manikins.

### **FEATURES**

- 1 | Soft and flexible genitalia allow proper cleaning and procedural preparation
- 2 | Successful insertion confirmation by fluid outflow
- 3 | Enema insertion in lateral position
- 4 | Detachable external genitalia

#### SKILLS

- Intermittent catheterization
- | Indwelling catheterization
- Perineal care
- I Enema
- | Localization of urethral meatus
- | Cleansing and aseptic technique
- | Catheter insertion and placement
- | Catheter Removal
- | Manual bladder compression

Compact sets consist of a genitalia unit and accessories, for hybrid training and table top demonstration.

11381-199-B

Male Genital Unit only

MW2D

11433-299-B

Female Genital Unit only







**Enema Training** 



Pass of catheter can be observed through transparent bladder reservoir



**DESCRIPTIONS** 

#### SPECIFICATIONS

Size: 34 x 33 x 24cm / 13.5 x 13 x 9.5 inch 
Soft resin 
Weight: 3.5kg / 7.7lbs(M) 
3.4kg / 7.5lbs(F) 
Latex free SET INCLUDES

- 1 lower torso manikin
- genital organ unit
- valves (one is incorporated in the assembled unit)
- 1 irrigation bottle

### MATERIALS

- rectum tube
- irrigation drainage tubes
- lubricant
- storage case
- instruction manual

#### REPLACEMENT PARTS

11381-100-01 clip (small) for urethral tube for MW2A 11381-100-02 clip (large) for MW2A/MW2B/MW2C

11433-000-01 clip (large) for MW52

11381-100-03 urinary bladder reservoir

11381-100-04 rectum unit

RECOMMENDED DEVICES

11381-300 6 replacement valves

11433-010 1 replacement female skin

11433-020 2 replacement female meatus

11381-800 lubricant

11381-410 1 replacement male skin

11381-510 2 replacement male meatus

\* Catheters are not included.