

Puerperal Uterus Palpation Training Model (Wearable type)

LM-116

This is a model that can be used to examine postpartum uterine involution, and can be worn by a person or a mannequin. The model represents the state of uterine involution up to 5 days postpartum. The model can be worn by a person to practice postpartum care, including communication training.



Features

- The abdominal wall is created from a special material to closely resemble the texture of a human body. It is well suited for demonstrating palpation, because the model allows the fundus to be examined by palpation.
- By changing the uterus models (four types), it can be used to check the status of uterine involution during the postpartum period.
- By having a person wear the model, postpartum care can be demonstrated while communicating with them.
- By wearing the model, a person can experience how a woman feels after giving birth.
- The model includes a realistic vulva to remind users that a postnatal mother might have feelings of modesty and embarrassment.
- The model can be used to demonstrate realistically how to change sanitary shorts and pads.

Specifications

Puerperal Uterus Palpation Training Model (Wearable type)

Compatible size

Waist 54cm to 116cm
(140cm with the extension pad)
Hip 54cm to 116cm
(140cm with the extension pad)

Total weight

Approx. 3.7kg

Storage bag

Length 31cm Width 46cm Height 22cm

Please note that product specifications are subject to change without notice to improve the product.

KOKEN CO., LTD.

INTERNATIONAL MARKETING & SALES
1-4-14 Koraku, Bunkyo-ku, Tokyo 112-0004, Japan
Tel: 81-3-3816-3562 © 2022 KOKEN CO., LTD.

www.kokenmpc.co.jp

