

Haemocytometer



GrayMed
PH: 61 3 86091732
Email: sales@graymed.com.au
Web: www.graymed.com.au

Our Haemocytometer is specially designed for counting erythrocytes and leukocytes. Double rulings are visible through a wide range of illumination. Neubauer rulings are 0.1mm below cover glass. The H-shaped moat offers two counting areas. Fused rhodium on glass provides distinct lines.



Haemocytometer

Model #	Description
CTL-HEMM-GLDR	Haemocytometer

FEATURES

- Bright-Line
- Double Neubauer Improved
- Designed for counting erythrocytes and leukocytes
- H-shaped moat offers two counting areas
- Neubauer rulings are 0.1 mm below cover glass

SPECIFICATIONS

NB. Manufactured in ISO 13485 certificate approved facilities according to the IVD 98/79/EC directive and authorised for supply in Australia by the TGA for use as an IVD medical device. See TGA regulatory framework for in-vitro diagnostic medical devices.

CONSUMABLES

Cover Slips. Product Code: CTP-CVR7-GL77
Additional cover slips for haemocytometer
Glass, 20x26mm, 10 Pieces per pack.



P 770.270.1394
F 770.270.2389

865 Marathon Parkway
Lawrenceville GA 30046

LWScientific.com