



DR. WACK
Made in Germany

Product-Information

The Problem:

Until now, there was no such thing as a care product for suede and matt, smooth leather that could both maintain their natural appearance and restore important nutrients in the leather. Without additional nutrients, suede and matt, smooth leather can become rough and brittle, and their colours can fade.

The Solution:



Product Advantages:

Permanent moisture protection and easy care for matt smooth leather and rough leather



- Keeps leather feeling soft and smooth
- Reduces dirt adhesion
- Water-repellent
- Colourless, but reinvigorates original colours
- Preserves the breathability of climate membranes

Applications:

Matt smooth leather (incl. TFL Cool, perforated leather) as well as all types of rough leather, such as nubuck, suede and velour leather

- Also suitable for leather/textile combinations



How to apply:

Shake well before use. Spray product evenly and not excessively onto the dry and clean fabric from a distance of approximately 20 cm. Leave to dry in a shaded place. Do not use in direct sunlight or on hot surfaces.

Tip:

Repeating the application process once the fabric has completely dried intensifies the level of protection (from moisture and dirt).



The product's liquid content must be two thirds full. Returned consignments must be correctly stamped and proof of purchase (invoice) included.