



Bell 407

Medical Interior

Phone: +1 (360) 225-1212
Email: solutions@lifeport.com

This document contains no technical data subject to the EAR or the ITAR

1610 Heritage Street
Woodland, WA. 98674

© 2020 LifePort



Bell 407 Medical Interior

Med Deck

Lightweight patient receiving deck pivots easily out of the aircraft door allowing for forward, aft and lateral movement.

Med Wall

Typical services include:

- 2 electrical outlets
- 1 air outlet with gauge and control switch
- 1 regulated vacuum outlet with gauge and control switch
- 1 oxygen outlet
- 1 regulated oxygen outlet
- 1 digital oxygen gauge
- Locking drug box
- Radio rack
- Netted storage
- Vacuum canister holder



Phone: +1 (360) 225-1212
Email: solutions@lifeport.com

This document contains no technical data subject to the EAR or the ITAR

MACHINED SLED

Lightweight machined aluminum stretcher with adjustable backrest and foot bag.

INCUBATOR CONVERSION

Compact with low-profile containment for oxygen and compressed air tanks (tanks sold separately) as well as customer-supplied units and IV-pumps.

MEDPAK

Kevlar oxygen supply, 1000-watt lightweight inverter, vacuum system, compressed air system, shore power outlet and door assembly with remote electrical shut-off.

PILOT BARRIER

Isolates the pilot compartment from the patient and protects the pilot from the AeroSled during loading.

BLACKOUT CURTAIN

Provides separation between the pilots compartment and medical team for NVG operations.

ACCESSORIES

- Device mounts: Several options for mounting monitors and other medical equipment.
- Portable Oxygen Bottle Mount: Lightweight mount used for extra oxygen cylinder. Composite fuel cell cabinet with radio console.

LIFEGARD FLOORING

Unique one-piece construction provides a seamless, sealed covering that is custom-formed to the aircraft floor.

1610 Heritage Street
Woodland, WA. 98674

© 2020 LifePort