


Features & Benefits

- ▶ Revolutionary curved design for ease of use
- ▶ Tapered blade front to ease the opening of ring-pull cans
- ▶ Perfectly balanced
- ▶ Spinner-ring finger protector
- ▶ Constructed from hand polished stainless-steel



SPECIFICATIONS	
Item Number	46/ENDOBLADE
Dimensions	
A - Length	185 mm (7.2 in)
B - Width	48 mm (1.88 in)
C - Height	14 mm (0.55 in)
D - Diameter	N/A
Weight	92 grms (0.20 lbs)
Construction	Stainless-Steel (hand polished)
Custom Logo Area	
E - Length	40 mm (1.57 in) (Logo must be positioned on the convex side)
Decoration Technique	Laser engraved
Important Note	Each EndoBlade™ is uniquely handcrafted and as such slight variations in finish will be apparent from one piece to another. This is an important aesthetic of this tool.

Des.Reg.Pat.

EndoBlade™ Anatomy	Usage, Washing & Caring
<p>1</p> <p>A. Spinner ring B. Bottle/ring-pull can opener</p>	<p>2</p> <p>Curve side facing up</p> <p>spinner ring</p> <p>Rear bottle opening edge</p> <p>To open bottles, place the Endoblades™ rear bottle opening edge under the crown seal of the bottle.</p> <p>Pull the Endoblade™ upwards to release the bottle top.</p> <p>Curve side facing up</p> <p>Pull the Endoblade™ upwards to open the can.</p> <p>Can tab opening. Place the tapered blade front under the ring pull tab.</p> <p>To facilitate quick draw, ensure that the Endoblade™ is positioned in your pocket with the curve side facing out.</p>
<p>WARNING NEVER CLEAN WITH ABRASIVE MATERIAL OR HARSH CLEANING CHEMICALS/SOLVENTS THOROUGHLY DRY AFTER CLEANING</p>	<p>WARNING NEVER PLACE ANY PART INTO A GLASS WASHER, DISH WASHER OR BOILING WATER</p> <p>The EndoBlade™ can be cleaned by wiping with a soft cloth that has been immersed in warm/soapy water</p>