



04

Cosmetics

Cosmetics

Home Bleaching	p.35
In-Office Bleaching	p.35
Non-Vital Walking Bleaching	p.36
White Spots Removers	p.36
Accessories	p.36 - p.37
Others	p.37

Cosmetics

Home Bleaching

Opalescence® PF 10%, 15%, 20%, and 35% / Ultradent

Carbamide Peroxide with Potassium Nitrate and Fluoride



- PF formula strengthens enamel, decreases sensitivity, prevents caries
- Formulated to prevent dehydration and shade relapse
- Sticky, viscous gel won't migrate to soft tissues and ensures tray stays securely in place
- Day or night wear



Suggested Opalescence® Wear Times:

- PF10%** - 8 to 10 hours or overnight
- PF15%** - 4 to 6 hours
- PF20%** - 2 to 4 hours
- PF35%** - 30 minutes

Home Bleaching

Opalescence® GO 10% Mini Kit / Ultradent ✓ Pre-Loaded Whitening Tray



**10%
Hydrogen Peroxide**



Wear 30 mins a Day Only!!!

The New UltraFit tray adapts quickly to any smile & provides molar-to-molar coverage



Before

After

Even Works Over Braces!



1 BOX for whitening touch-up

2 BOXES for a standard course of whitening

3 BOXES for severe discoloration

Flavor	Packing	Item no.
Mint	4 pairs/box	UL-4645-2

■ Patient Kits

- Contains: - 8 x Syringes (1.2 mL)
- 1 x Opalescence whitening toothpaste (1 oz.)
- 1 x Tray case
- 1 x Shade guide and carry all bag

Flavor	10%	10%PF	15%PF	20%PF	35%PF
Mint	UL-5361	UL-5364	UL-5367	UL-5370	UL-5373
Melon	UL-5362	UL-5365	UL-5368	UL-5371	UL-5374
Regular	UL-5363	UL-5366	UL-5369	UL-5372	UL-5375

■ Doctor Kits

- Contains: - 8 x Syringes (1.2 mL)
- 1 x LC Block-Out Resin syringe (1.2 mL) and Black Mini tip
- 2 x Sof-Tray 0.035" (5" x 5") sheets
- 1 x Opalescence whitening toothpaste (1 oz.)
- 1 x Tray case
- 1 x Shade guide and carry all bag

Flavor	10%	10%PF	15%PF	20%PF	35%PF
Mint	UL-5376	UL-5379	UL-5382	UL-5385	UL-5388
Melon	UL-5377	UL-5380	UL-5383	UL-5386	UL-5389
Regular	UL-5378	UL-5381	UL-5384	UL-5387	UL-5390

■ Refills

- Contains: - 40 x Syringes (1.2 mL)
- 1 x Shade guide

Flavor	10%	10%PF	15%PF	20%PF	35%PF
Mint	UL-5391	UL-5394	UL-5397	UL-5400	UL-5403
Melon	UL-5392	UL-5395	UL-5398	UL-5401	UL-5404
Regular	UL-5393	UL-5396	UL-5399	UL-5402	UL-5405

In-Office Bleaching

Opalescence® BOOST PF 40% / Ultradent



- **No light needed!**
- **No equipment use!**
- Syringe-syringe mix - ensures freshness
- increase micro hardness of enamel, reduce caries, minimize sensitivity
- **Use:**
20 minutes per session
2-3 sessions per treatment

40% Hydrogen Peroxide with Potassium Nitrate & Fluoride



Before

After

■ Introductory Kit

Content	Item no.
- 4 x Opalescence Boost / Activator syringes (1.2 mL) - 2 x OpalDam Green syringes (1.2 mL) - 2 x Ultradent Luer Vacuum Adapters	UL-4750
- 2 x IsoBlocks	
- 10 x Micro 20ga tips	
- 2 x Shade guide card - 10 x Micro 20ga FX tips	
- 2 x SurgiTips	

■ Patient Kit

Content	Item no.
- 2 x Opalescence Boost / Activator syringes (1.2 mL) - 1 x OpalDam Green syringes (1.2 mL) - 1 x IsoBlocks	UL-4751
- 5 x Micro 20ga tips	
- 1 x Shade guide card - 5 x Micro 20ga FX tips	

■ Refill

Content	Item no.
- 2 x Opalescence Boost / Activator syringes (1.2 mL) - 5 x Micro 20ga tips	UL-4753
- 1 x Shade guide card	



Non-Vital Walking Bleaching

Opalescence® Endo / Ultradent

- Formulated to whiten endodontically treated, discolored teeth
- Easy to place inside pulp chamber
- 3-5 day treatment

35% Hydrogen Peroxide



Mini Kit

Content	Item no.
- 2 x Syringes (1.2 mL) - 20 x Black Mini tips	UL-1270

Mini Refill

Content	Item no.
- 2 x Syringes (1.2 mL)	UL-1323

White Spots Removers

Icon® / DMG

✓ Caries Infiltration

- Esthetic results on smooth surfaces
- Preservation of heating tooth structure
- Pain-free method without anesthesia or drilling

For Proximal Treatment - Intro Kit



Content	Item no.
2 x Treatment units each including: - 1 x Syringe @ 0.3 ml Icon-Etch - 1 x Syringe @ 0.45 ml Icon-Dry - 1 x Syringe @ 0.45 ml Icon-Infiltrant - 6 x Approximal-Tips - 4 x Dental wedges - 1 x Luer-Lock-Tip	DMG-220341

For Smooth Surface Vestibular Treatment - Intro Kit



Content	Item no.
2 x Treatment units each including: - 1 x Syringe @ 0.3 ml Icon-Etch - 1 x Syringe @ 0.45 ml Icon-Dry - 1 x Syringe @ 0.45 ml Icon-Infiltrant - 6 x Vestibular-Tips - 1 x Luer-Lock-Tip	DMG-220343

Opalustre® / Ultradent

✓ Chemical and Mechanical Abrasion Slurry

- Quickly removes superficial white / brown stains 0.2mm in depth

6.6% Hydrochloric Acid Slurry



Kit

Content	Item no.
- 4 x Syringes (1.2 mL) - 20 x White Mac Tips - 20 x OpalCups (10 Bristle & 10 Finishing)	UL-554


Refill

Content	Item no.
- 4 x Syringes (1.2 mL)	UL-555

Accessories

Tray Sheets

HENRY SCHEIN®



Vacuum-Forming Materials

- Thermoforming Sheets
- Size: 5" x 5"

Base Plate Material

Description	Color	Thickness	Packing (Sheets per Box)	Item no.
Hard	Pink	0.060"	25 sheets/box	HS-100-4623
			625 sheets/box	HS-100-2306

Coping Material

Description	Color	Thickness	Packing (Sheets per Box)	Item no.
Soft	White	0.020"	50 sheets/box	HS-100-3513
		0.030"		HS-101-5448

Temporary Splint Material

Color	Thickness	Packing (Sheets per Box)	Item no.
Clear	0.020"	50 sheets/box	HS-100-3822

Mouth Guard Material

Color	Thickness	Packing (Sheets per Box)	Item no.
Clear	0.040"	12 sheets/box	HS-101-8910
	0.060"		HS-101-7962
	0.080"		HS-100-7718
	0.120"		HS-101-2123
	0.150"		HS-100-1199
Red	0.150"	12 sheets/box	HS-100-7408
Blue			HS-101-8361
White			HS-101-8270
Yellow			HS-101-4479

Surgical Tray Material

Color	Thickness	Packing (Sheets per Box)	Item no.
Clear	0.060"	25 sheets/box	HS-100-1330
	0.080"		HS-100-4239
	0.100"		HS-101-4413

Tray Material

Color	Thickness	Packing (Sheets per Box)	Item no.
Pink	0.080"	25 sheets/box	HS-100-2184
	0.100"		HS-100-6962
	0.125"		HS-100-2195
Blue	0.150"		HS-100-6997



Anaesthetics
Articulating
Bur and Abrasives
Cosmetics
Diagnostic
Disposables
Endodontics
Equipment
Handpieces
Implant
Infection Control
Instruments
Isolation
Laboratory
Oral Surgery
Orthodontics
Patient Education and Gift
Pedodontics
Periodontal
Post and Core
Preventive
Prosthodontics
Restorative
Uniform

Cosmetics

Accessories

Tray Sheets

Sof-Tray® Classic Sheets

/ Ultradent



Size	Thickness	Packing	Item no.
Sof-Tray Regular Sheets (For most bleaching trays)			
5"x 5" (127mm x 127mm)	0.035" (0.9mm)	25 sheets/box	UL-226
Sof-Tray Medium Sheets (For bleaching patients who are bruxers)			
5"x 5" (127mm x 127mm)	0.060" (1.5mm)	20 sheets/box	UL-227
Sof-Tray Heavy Sheets (For severe bruxers or a TMJ appliance)			
5"x 5" (127mm x 127mm)	0.080" (2.0mm)	20 sheets/box	UL-284

Others

Liquid Dam

OpalDam® and OpalDam® Green / Ultradent

- Protects soft tissue with impervious seal
- Removes easily
- Applies directly

Light-Cured Resin Barriers



OpalDam - A passively adhesive (sealing) methacrylate-based resin barrier used for isolating tissue adjacent to teeth being whitened

OpalDam Green - Ensures a safe, unmistakable barrier every time

Kit

Content	Item no.
- 4 x Syringes (1.2 mL) - 10 x Black Mini Tips	UL-324

Kit

Content	Item no.
- 4 x Syringes (1.2 mL) - 10 x Black Mini Tips	UL-1824

Refill

Content	Item no.
4 x Syringes (1.2 mL)	UL-325

Refill

Content	Item no.
4 x Syringes (1.2 mL)	UL-1825

Econo Refill

Content	Item no.
20 x Syringes (1.2 mL)	UL-326

Econo Refill

Content	Item no.
20 x Syringes (1.2 mL)	UL-1826

Block Out Resin

LC Block-Out Resin / Ultradent

Low Viscosity



- Provides reservoir space for bleaching trays
- Blue pigment and great utility resin

Kit

Content	Item no.
- 4 x Syringes (1.2 mL) - 20 x Black Mini Tips	UL-240

Refill

Content	Item no.
4 x Syringes (1.2 mL)	UL-241

Econo Kit

Content	Item no.
- 20 x Syringes (1.2 mL) - 20 x Black Mini Tips	UL-242

Tray Trimming Scissors

Utility Vinyl Cutters / Ultradent

- Use for gross trimming of tray

Packing	Item no.
Each	UL-604



Ultra-Trim Scalloping Scissors / Ultradent

- Use for precision trimming of border around interdental papilla
- Stainless steel



Packing	Item no.
Each	UL-605





