

# Wireless Algometer

COMMANDER<sup>®</sup>  
ECHO<sup>™</sup>

**Wireless Commander Echo Algometer easily measures pressure thresholds and pain tolerances reported by patients. Algometry is ideal for helping evaluate fibromyalgia and trigger point tenderness.**



## Advanced Features

- ▶ Innovative pop-in “quick connect” system makes changing test tips fast and easy
- ▶ Console stores 20 tests with up to 4 repetitions per side
- ▶ Quantify joint tenderness in arthritic conditions
- ▶ Evaluate pain sensitivity and identify trigger points
- ▶ Thin profile improves stability during tests
- ▶ Ergonomic shape fits your hand and provides a secure grip
- ▶ Wireless console is independent of the testing device allowing the tester to review data during active testing without influencing the patient’s level of effort



\*Console required to operate device. Console and device sold separately. Device shown with accessory attached.

DISTRIBUTED BY



190, Boul. St-Elzéar Ouest  
Laval (Québec) H7L 3N3  
T 450 669.8985 | 1 800 361.0877  
F 450 669.9532 | 1 888 935.7001  
[www.medicaltronik.ca](http://www.medicaltronik.ca)



**JTECH MEDICAL**  
MEASURING HUMAN PERFORMANCE  
© 2012 JTECH MEDICAL MM320\_C

# Technical Specifications

Number of tests stored on Console: . . . . . 20  
 Number of repetitions per test: . . . . . 1-4  
 Audible cues: . . . . . yes  
 Load cell capacity: . . . . . 25 lb.  
 Auto coefficient of variation calculations (CV):. yes  
 Automatic bilateral deficit calculation:. . . . . yes  
 Threshold settings: . . . . . 0.2, 0.4, 1, 3 lb.  
 Units of measure (pounds default): . . . lb, kg, N  
 Automatic fatigue calculation:. . . . . yes  
 Timed test for fatigue evaluation . . . . . yes  
 Allows for specific test re-testing:. . . Last rep only while still testing

Print options: . . . . . Yes\*  
 \* With purchase of Downloader Echo software

### Included Accessories:


**0.5cm<sup>2</sup> Algometer Tip** (8AC006) – Can be used for cervical testing.

**1.0cm<sup>2</sup> Algometer Tip** (8AC005) – Conforms to normative testing.

**Calibration Flat Pad** (8AC008) – Used to assist in calibrating device.

### Optional Accessories:

**Finger Tip Adapter** (8AC007) – Converts the Algometer for testing individual finger strength.



## Product Highlights

- Console records and stores up to 20 tests
- Fine resolution that identifies small yet clinically significant changes in pressure sensitivity
- Replaceable rubber tips
- Two tips (1cm<sup>2</sup> and 0.5cm<sup>2</sup>) for applying to either large muscles and small muscles of the head and neck
- Document long and short term effectiveness of treatment
- “Smart” load cell technology maintains accuracy even when force is applied off-center.

## COMMANDER ECHO Downloader™

Optional Commander Echo Downloader (CM311) software enables you to create quick but meaningful reports from the data you collected on your Commander Echo Console with the click of a button. The easy-to-use Commander Echo Downloader software interfaces directly with your Commander Echo Console eliminating dictation and transcription errors.



Other  
Commander  
Echo Devices\*



\*Some devices shown with accessories attached.

DISTRIBUTED BY



190, Boul. St-Elzéar Ouest  
 Laval (Québec) H7L 3N3  
 T 450 669.8985 | 1 800 361.0877  
 F 450 669.9532 | 1 888 935.7001  
[www.medicaltronik.ca](http://www.medicaltronik.ca)



**JTECH MEDICAL**  
 MEASURING HUMAN PERFORMANCE

© 2012 JTECH MEDICAL MM320\_B