

zeta 3 wipes TOTAL

Alcohol solution wipes for the rapid disinfection and cleaning of the surfaces of small medical devices.



Zeta 3 total wipes Refill (120 wipes bag refill)
C810063 Zeta 3 Wipes Dispenser (120 wipes tub)

zeta 3 wipes POP-UP

Large wipes soaked in an alcohol-free solution for the rapid disinfection and cleaning of medical device surfaces, including delicate surfaces (dental units in real or artificial leather, plexiglass surfaces, etc.).



C810064 Zeta 3 Wipes POP-UP
100 wipes soft pack

zeta 6 hydra

An ultra delicate detergent for frequent use on hands and sensitive skin. Gentle on the skin and enriched with skinprotecting substances and nutrients. Zeta 6 Hydra does not contain SLES and Parabens found in conventional soap.



C810042 Zeta 6 Hydra 1 litre bottle with dosage cap
C810043 Zeta 6 Hydra 5 litres container with dosage cap

zeta 1 ultra

A concentrated, broad spectrum, liquid disinfectant and detergent. Suitable for all surgical and rotary instruments, even the most delicate (scalpels, tweezers, burs, mirrors, probes, etc.).



C810000 Zeta 1 Ultra 1 litre bottle
C810001 Zeta 1 Ultra 5 litres container with dosage cap

zeta 3 foam

A ready-to-use, alcohol-free disinfectant and detergent foam. Especially suitable for more delicate medical device surfaces (handpieces, contra-angles, dental units, etc.).



C810025 Zeta 3 Foam
750 ml bottle with nozzle trigger spray

C810026 Zeta 3 Foam
3 litres container with dosage cap


zeta 7 spray & solution

Zhermack offers excellent performance products for impression disinfection. A complete spectrum (bactericidal, fungicidal, tuberculocidal and virucidal) and excellent compatibility with different types of impression material.



1% DILUTION

C810048 Zeta 7 Solution 1 litre bottle
C810150 Zeta 7 Spray 750 ml



zeta 5 power act

The 3 in 1 solution to disinfect, clean and deodorise aspiration circuits



C810040 1 L bottle

1% DILUTION

One litre of ZETA 5 Power Act gives 100 litres of disinfectant solution

▶ EFFECTIVE

Thanks to its broad spectrum of action, Zhermack disinfectant and detergent guarantees high protection and safety for both user and patients.

▶ PRACTICAL

It satisfies the surgery's needs and schedules: as well as for standard evening disinfection*, it can also be used for a rapid** 15 minute action between patients, or after particularly invasive operations.

▶ SAFE

For the user, as it is aldehyde and phenol free. For the aspiration circuits, thanks to its non-foaming formula which helps keep them in good condition and guarantees regular functioning.

*bactericidal, yeasticidal, tuberculocidal, limited virucidal (for example against HIV, HBV and HCV) ** excluding TBC action



Bio Medical Surgical Equipment & Instruments Trading L.L.C.
Tel: +971 4 3596242,
P.O. Box: 12134, Dubai U.A.E.
Website: www.biolinkdubai.com
E-mail: info@biolinkdubai.com



Zhermack S.p.A. Via Bovazecchino, 100
45021 Badia Polesine (RO) ITALY
T +39 0425 597611 | F +39 0425 53596
info@zhermack.com | www.zhermack.com

elite glass

- High transparency of the material
- Excellent photopolymerization of the composite used for restoration
- High quality standards of the aesthetic restoration
- Versatility of use



C401610 Standard Pack:
2 x 50 ml cartridges + 6 mixing tips green

Products information reference

DUBAI



Algitray



Zeta 1 Ultra



Zeta 2 Sporex

hydrorise

High performance for Hydrorise, the hyperhydrocompatible A-silicone for high precision impression taking, which adapts to different impression techniques (double step impression, simultaneous impression with dual viscosity and monophasic impression) in various clinical situations.




Hydrorise Putty
C207011 **Standard Pack:** 300 ml Base + 300 ml Catalyst tubes
C207013 **Eco Pack:** 900 ml Base + 900 ml catalyst tubes



Hydrorise Maxi Heavy Body
C207043 2 x 380 ml cartridges + 15 dynamic-static mixing tips + 2 mixing tip lockers
Hydrorise Heavy Body
C207009 **Standard Pack:** 2 x 50 ml cartridges (Base+Catalyst) + 6 mixing tips green
Hydrorise Light Body
C207001 **Standard Pack:** 2 x 50 ml cartridges + 12 mixing tips yellow

Product	Delivery system	Viscosity	Recommended techniques	Implantology	Type of setting	Working time* (min:s)	Time in mouth* (min:s)	Setting time* (min:s)
Hydrorise Implant Heavy Body	5:1 automatic mixing	Heavy	Single step impression	●	Normal set	02:00	03:30	05:30
Hydrorise Implant Medium Body	5:1 automatic mixing	Medium	Monophasic impression	●	Normal set	02:00	03:30	05:30
Hydrorise Implant Light Body	1:1 dispenser	Light	Single step impression	●	Normal set	02:00	03:30	05:30

* the times mentioned must be intended at 35 °C - 95 °F



VPS Tray Adhesive



AcryTemp



Hydrorise Light Body
C207001 **Standard Pack:** 2 x 50 ml cartridges + 12 mixing tips yellow



Hydrorise Maxi Monophasic
C207041 2 x 380 ml cartridges + 15 dynamic-static mixing tips + 2 mixing tip lockers



Hydrorise Monophasic
C207007 **Standard Pack:** 2 x 50 ml cartridges (Base+Catalyst) + 6 mixing tips

zetaplus system

High performance condensation silicones for impressions in fixed and removable prostheses.

zetaplus

The «flow conditioning» effect helps the material flow and reduces pressure on the soft tissue. Due to its rigid set, Zetaplus is recommended putty for the two-stage technique.



C100600 **Standard pack:** 1 tube 900 ml (1.53 kg)



indurent gel

Gel catalyst for C-Silicones
In combination with Zetaplus System



C100700
1 tube 60 ml

oranwash L

This light bodied hydrophilic flowable silicone is suitable for the two-stage technique.



C100660 **Standard pack:**
1 tube 140 ml

hydrogum 5

- High dimensional stability
- Quick working and setting times
- Thixotropy
- Pleasant scent
- Long dimensional stability of impressions (5 days of dimensional stability)
- Easy mixing
- Homogeneous mass, smooth, compact surfaces

C302070
Hydrogum 5 Extra Fast Set 1 bag 453 g



tropicalgin

- Fast working and setting times
- Thixotropy
- Chromatic alginate which change colour to guide the dentist during the 3 processing stages
- Easy mixing
- Homogeneous mass, smooth, compact surfaces

C302240
Tropicalgin Fast Set 1 bag 453 g



orthoprint

- Highly elastic alginate
- Extra-fast working and setting times
- Vanilla scent
- Fast water absorption
- Easy mixing
- Emetic reflex reduction

C302145
Orthoprint Extra Fast Set 1 bag 500 g



occlufast rock

- A-silicones for bite registration
- Imperceptible consistency
- Minimum permanence in the oral cavity
- Thixotropic: allows you to easily control the positioning on the occlusal plane
- High final hardness (95 Shore A)

C200726
Occlufast Rock, 2 x 50 ml cartridges + 12 mixing tips green



colorbite D

- Thermochromic: a technology which uses temperature sensitive pigments
- Lime scent, green color at 23 °C (73 °F) and yellow over 30 °C (86 °F)
- Quick setting in the oral cavity
- Extremely rigid and hard (45 Shore D)

C200750
Colorbite D, 2 x 50 ml cartridges + 12 mixing tips green

