# As used by the NHS Deep Cleaning and Advisory Service.

Derived from the electrolysis of organic salt and water and it is the preferred cleaning agent for hard surfaces and soiled upholstery before disinfection with Pathisol Clinical Cleaner and Disinfectant.

#### The benefits are clear

- + Use as supplied, no dilution needed
- A universal cleaning solution that lifts dirt from virtually any surface and is as gentle on the environment as water
- Used in live environments
- Can be used in existing applicators/machinery
- Helps fund frontline patient care for NHS Trusts
- Available in 600ml to 1000ltr containers
- Easy to apply and works immediately. Sometimes the more determined dirt needs repeated applications and vigorous wiping
- + Delivered to your door each month
- Application training provided in line with NHS standards

#### What makes Pathisol unique?

## Safe to use around people and works in seconds

Pathisol Heavy Duty All Purpose Cleaner is an effective degreasing and cleaning solution. Apply liberally to the area to be cleaned and wipe up with a microfibre cloth. We simply deliver a fresh supply of the required volume to your door each month.

# Directly supports frontline patient care for the NHS

Pathisol Heavy Duty All Purpose Cleaner is a hospital made degreasing and cleaning solution produced by Tutum Health at the University Hospital of Hartlepool. The income generated by Pathisol enables reinvestment into NHS Trusts' to enhance frontline patient care and services.

#### For enquiries and advice, please call 0808 164 2800





# pathisol all purpose cleaner and degreaser

Hospital made from organic ingredients

## Additional infection prevention solutions and services

Working with schools, businesses and other organisations to bring infection prevention practices in line with NHS standards:

- + Site auditing and retraining of existing cleaning staff with official certification
- + Deep cleaning and biodecontamination with official certification all starts with Pathisol All Purpose Cleaner being used as the cleaning solution before the application of the Pathisol disinfectant for the best results
- + Sustainable protective equipment solutions

