

www.antibodiesinc.com orders@antibodiesinc.com 530-758-4400

## **Product Datasheet**

## Anti-GABA-A-R, Alpha1 Antibody FL650 Conjugate



## Overview

Catalog #	75-136-FL650
Conjugate	FL650 Ex: 655 nm, Em: 676 nm
Isotype	lgG2a
Clone Number	N95/35
Size	200 μL
Concentration	0.5 mg/mL
Host Species	Mouse Monoclonal
Format	Purified by Protein A chromatography
Buffer	PBS with 0.09% azide
Applications	ICC, IHC
Species Reactivity	Human, Mouse, Rat, and Xenopus
Immunogen	Fusion protein amino acids 355-394 of mouse GABA-AR-Alpha1 (accession number P62812) produced recombinantly in E. Coli
Molecular Weight	55 kDa
Cite this Antibody	Antibodies Inc Cat# 75-136-FL650, RRID: AB_2939507
Details	
Details Target Description	GABA(A)R, Alpha1 (Gamma-aminobutyric acid receptor subunit alpha-1) makes up two of the five subunits in the membrane bound GABA A receptor, the others being two β2 subunits and a single γ2 subunit. The receptor is a ligand gated chloride ion channel found postsynaptically and binds GABA, the major inhibitory neurotransmitter in the CNS. Mutations in GABA(A)R, Alpha1 have been linked to several types of epilepsy.
	subunits in the membrane bound GABA A receptor, the others being two $\beta$ 2 subunits and a single $\gamma$ 2 subunit. The receptor is a ligand gated chloride ion channel found postsynaptically and binds GABA, the major inhibitory neurotransmitter in the CNS. Mutations in GABA(A)R, Alpha1 have been
Target Description	subunits in the membrane bound GABA A receptor, the others being two $\beta$ 2 subunits and a single $\gamma$ 2 subunit. The receptor is a ligand gated chloride ion channel found postsynaptically and binds GABA, the major inhibitory neurotransmitter in the CNS. Mutations in GABA(A)R, Alpha1 have been linked to several types of epilepsy.
Target Description Specificity	<ul> <li>subunits in the membrane bound GABA A receptor, the others being two β2 subunits and a single γ2 subunit. The receptor is a ligand gated chloride ion channel found postsynaptically and binds GABA, the major inhibitory neurotransmitter in the CNS. Mutations in GABA(A)R, Alpha1 have been linked to several types of epilepsy.</li> <li>No cross-reactivity against GABA-A-R-Alpha2 or -Alpha3</li> <li>Produced by in vitro bioreactor culture of hybridoma line followed by Protein A affinity</li> </ul>

## **Our Guarantee**

As an original manufacturer, we are dedicated to creating quality and reproducible antibodies that further your research. We provide personalized customer support from the scientists that made the antibody and offer a free replacement or 100% refund if we cannot resolve an issue. Order today and experience our 50+ year passion for science.

Note: For research use only. Not intended for therapeutic or diagnostic use. Use of all products is subject to our terms and conditions, viewable on our website.