

Dermastir Ampoules - Luxury Serum

Dermastir Ampoules are luxurious anti-ageing serums that contain concentrated active ingredients. The active ingredients provide targeted action against signs of ageing and various skin conditions, improving the quality of the skin and its aspect.

Dermastir Ampoules are presented in single-dose 2-3 ml ampoules made of special medical glass that prevents any risk of oxidation or contamination. Thanks to this monodose packaging, the serums deliver to the skin the freshest formula with maximum performance of active ingredients.


SERUMS

MONODOSE

Medical glass ampoules



Quantity: 10 x 3ml



Dermastir Ampoules can be used in combination with aesthetic machines such as ultrasound, galvanic, high frequency, RF, PDT LED

The high concentration of active ingredients of Dermastir Ampoules provides visible results which last even longer if ampoules are used as course treatments. Dermastir Ampoules can be used at home, during aesthetic treatments and with aesthetic machines. Dermastir Ampoules can be used directly on the skin or mixed with a cream or the mask.

DERMASTIR®

WWW.DERMASTIR.COM



LABORATOIRES
ALTA
CARE

6, Piazza di Spagna - 00187 Roma, Italia.
www.altacare.com

Dermastir Luxury Ampoules Collagen

The Dermastir Skincare Ampoules with Collagen are essential to improve skin elasticity and firmness. Collagen is a natural protein that provides structural support. It is found throughout the body in the skin, muscles, tendons and bones. Fibers of collagen are woven together like threads in fabric to form a framework into which new cells can grow. In the skin, collagen provides texture, resiliency, and shape. Human beings are constantly producing collagen, but as we age, our collagen breaks down faster than it can be replaced.

Dermastir Ampoules with Collagen improve the circulation and the tone of the skin and reinforce its elasticity and firmness.



FACIAL LIFT



ANTI-WRINKLE



DULL SKINTONE



UNEVEN SKINTONE



DIRECTION FOR USE



Break the top of the skincare ampoule. Distribute the precious serum directly on the skin or mix in a cream or in a mask. Massage delicately on the interested areas until completely absorbed.

NORMAL PROGRAMME

1. Apply one skincare ampoule every seven days for a month.
2. To consolidate the obtained results repeat phase 1.

INTENSIVE PROGRAMME


1. Apply one skincare ampoule every three days for a month.
2. To consolidate the obtained results repeat phase 1.

WARNINGS: Do not expose yourself to the sun. Avoid contact with eyes and in case rinse with water. The skincare ampoule tip must be broken by using a piece of cloth.

Collagen


One of the most well-known benefits of collagen is its ability to promote glowing, vibrant skin.

This essential protein provides elasticity to the skin, helping it to appear more youthful and healthier. But as we get older, collagen production declines, and fine lines, sagging skin, and dryness can occur.



Younger skin

High Density of Collagen and Elastin



Mature skin

Lack of Collagen and Elastin leads to loss of firmness

