

Dermastir Ampoules - Luxury Serum

Dermastir Ampoules are luxurious anti-ageing serums that contain concentrated active ingredients. The active ingredients provide targeted action against signs of ageing and various skin conditions, improving the quality of the skin and its aspect.

Dermastir Ampoules are presented in single-dose 2-3 ml ampoules made of special medical glass that prevents any risk of oxidation or contamination. Thanks to this monodose packaging, the serums deliver to the skin the freshest formula with maximum performance of active ingredients.


SERUMS

MONODOSE

Medical glass ampoules



Quantity: 10 x 3ml



Dermastir Ampoules can be used in combination with aesthetic machines such as ultrasound, galvanic, high frequency, RF, PDT LED

The high concentration of active ingredients of Dermastir Ampoules provides visible results which last even longer if ampoules are used as course treatments. Dermastir Ampoules can be used at home, during aesthetic treatments and with aesthetic machines. Dermastir Ampoules can be used directly on the skin or mixed with a cream or the mask.

DERMASTIR®

WWW.DERMASTIR.COM



LABORATOIRES
ALTA
CARE

6, Piazza di Spagna - 00187 Roma, Italia.
www.altacare.com

Dermastir Luxury Ampoules Glycolic Acid

Glycolic Acid is a naturally occurring acid that makes part of the Alpha Hydroxy Acids (AHAs). Dermastir Ampoules Glycolic Acid have delicate and at the same time potent action to dissolve the bonds between dead cells and provide gentle exfoliation to the skin surface. Dermastir Glycolic Acid Ampoules deal with pigmentation, dark spots, and signs of aging, reducing wrinkles and fine lines.

Dermastir Ampoules Glycolic Acid are indicated for oily skin. These ampoules are ideal both for young skin to restore brightness and freshness during cold months and tone the skin; and for mature skin that needs a fresh wave.



ANTI-WRINKLE



PIGMENTATION



DULL SKINTONE



Glycolic Acid

Exfoliating action favours the elimination of dead cells from the skin surface and is helpful in case of cutaneous dyschromia caused by sun-exposure and photoaging.

Beneficial effects on the skin prone to blackheads, white spots, blemishes and other imperfections.

Important astringent properties play role to diminish enlarged and blocked pores. Contribution to more even skin tone in case of xerosis, ichthyosis and wrinkles.

Glycolic acid stimulates collagen synthesis and helps to improve cell proliferation for a more youthful skin aspect.

The molecule size facilitates an easy and rapid absorption to the skin, providing the optimal results.

DIRECTION FOR USE



Break the top of the skincare ampoule. Distribute the precious serum directly on the skin or mix in a cream or in a mask. Massage delicately on the interested areas until completely absorbed.

NORMAL PROGRAMME

1. Apply one skincare ampoule every seven days for a month.
2. To consolidate the obtained results repeat phase 1.

INTENSIVE PROGRAMME

1. Apply one skincare ampoule every three days for a month.
2. To consolidate the obtained results repeat phase 1.

WARNINGS: Do not expose yourself to the sun. Avoid contact with eyes and in case rinse with water. The skincare ampoule tip must be broken by using a piece of cloth.

