

Dermastir Ampoules - Luxury Serum

Dermastir Ampoules are luxurious anti-ageing serums that contain concentrated active ingredients. The active ingredients provide targeted action against signs of ageing and various skin conditions, improving the quality of the skin and its aspect.


Dermastir Ampoules are presented in single-dose 2-3 ml ampoules made of special medical glass that prevents any risk of oxidation or contamination. Thanks to this monodose packaging, the serums deliver to the skin the freshest formula with maximum performance of active ingredients.

SERUMS MONODOSE

Medical glass ampoules



Quantity: 10 x 2ml



Dermastir Ampoules can be used in combination with aesthetic machines such as ultrasound, galvanic, high frequency, RF, PDT LED

The high concentration of active ingredients of Dermastir Ampoules provides visible results which last even longer if ampoules are used as course treatments. Dermastir Ampoules can be used at home, during aesthetic treatments and with aesthetic machines. Dermastir Ampoules can be used directly on the skin or mixed with a cream or the mask.

DERMASTIR®

WWW.DERMASTIR.COM



6, Piazza di Spagna - 00187 Roma, Italia.
www.altacare.com

Dermastir Luxury Ampoules Hyaluronic

Dermastir Luxury Skincare Ampoules Hyaluronic help to provide high levels of hydration and increase the volume of the skin. The ampoules are recommended especially for the sensitive and dry skin. Hyaluronan (also called hyaluronic acid or hyaluronate) is a non-sulfated glycosaminoglycan distributed widely throughout connective, epithelial, and neural tissues. It is one of the chief components of the extracellular matrix, and it contributes significantly to cell proliferation and migration. The average 70-kg man has roughly 15 grams of hyaluronan in his body, one-third of which is turned over (degraded and synthesized) every day. Hyaluronan is also a major component of skin where it is involved in tissue repair.



ANTI-WRINKLE



DULL SKINTONE



DRY SKIN



Sodium Hyaluronate

Since hyaluronic acid is naturally found in the dermis of the skin, products that contain sodium hyaluronate help to replace the hyaluronic acid that is lost as a result of the natural skin aging process and sun exposure. This results in the depletion of the epidermis, thinning, aging and decreased moisture in the skin. When applied topically Dermastir Hyaluronic will penetrate the skin's surface and replenish the natural levels of hyaluronic acid already present. It holds 1,000 times its weight in water. That is why it is such an effective moisturizer: when applied to the skin, it helps skin retain moisture, helps with wound healing, promotes better circulation and nutrient absorption. Results in firmer, plumper skin, which makes wrinkles less visible.

DIRECTION FOR USE



Break the top of the skincare ampoule. Distribute the precious serum directly on the skin or mix in a cream or in a mask. Massage delicately on the interested areas until completely absorbed.

NORMAL PROGRAMME

1. Apply one skincare ampoule every seven days for a month.
2. To consolidate the obtained results repeat phase 1.

INTENSIVE PROGRAMME

1. Apply one skincare ampoule every three days for a month.
2. To consolidate the obtained results repeat phase 1.

WARNINGS: Do not expose yourself to the sun. Avoid contact with eyes and in case rinse with water. The skincare ampoule tip must be broken by using a piece of cloth.

