


TIMS MVP™ 4.0

The ultimate solution for recording, editing, & archiving of MBS, FEES, & endoscopy studies



**Still using just fluoro or low quality DVD-video?
Upgrade to high resolution, improved efficiency, & synced audio.**

TIMS MVP™ is the fourth generation medical imaging product from TIMS Medical. Built for speed and productivity, it is designed to streamline your workflow and result in more studies per day from your imaging suites. Just some of the advanced features include:

High Resolution Video with Synced Audio

Automated recording & 4 times more detail than DVD-video

Label Video Clips During Recording

Save at least 5 minutes per study!

Send Studies in DICOM to PACS/VNA

Easy archiving & PACS disk space savings of 75% to 90%

Remote Review & Analysis with TIMS Review™


Frees up imaging suite for more studies per day

- Windows 10 Pro 64 bit workstation
- 30 frames per second
- Select patients from your worklist for automated input
- Triggers with your modality for automated recording
- Unlimited video recording length
- Your work is automatically saved at all times
- Custom label lists
- Archive in DICOM to any PACS/VNA
- Archive to CD/DVD, USB, & network
- Eliminate removable media from workflow – eliminate HIPAA risk

- Export studies in MP4
- Remote review & analysis with TIMS Review
- Graphic telestration & annotation
- Stopwatch & study timers
- Non-destructive video editing
- Comprehensive Support & Maintenance
- Built-in & online training documents
- User interface built for speed & productivity

Go to www.tims.com for more details & video demonstrations!

TIMS MVP & TIMS Review



Workflow: Speed & Productivity!

With over 4000 installations worldwide, TIMS Medical has continuously worked with SLPs and endoscopists to optimize their workflows. TIMS MVP 4.0 takes productivity to the next level for the power user, the new user, and everyone in between.

For the simplest of workflows:

- 1 - Input the patient information by selecting from the patient list in the worklist.
- 2 - Record the entire study with synced audio automatically.
- 3 - Quickly label your video series on-the-fly to save at least 5 minutes on every MBS study.
- 4 - Archive in DICOM to PACS/VNA and/or send to a TIMS Review system in your lab or office for analysis.

Extensive tools are available for quick & simple editing, annotating or telestrating on video, timing analysis, video review & analysis, and exporting in MP4 format for presentations and reports.

TIMS Review: TIMS Review is a software only version of the TIMS MVP product (no video acquisition hardware). One license of TIMS Review is included with every TIMS system package. TIMS Review can be installed on any Windows 10 computer in the speech pathology or endoscopy office/lab area. To speed workflow and immediately free up the fluoroscopy or endoscopy suite for the next patient, the study can be sent immediately over the hospital network from the procedure room to one or more TIMS Review systems throughout the hospital. Clinicians can then edit, annotate, review, and analyze the studies from their offices, while freeing up the procedure rooms to enable more studies per day for the hospital.

TIMS MVP Ordering Guide

TIMS 2000 EN Package (no audio)

Includes 24" LCD monitor, standard onsite installation & training, Video Isolator, Trigger Kit, & TIMS Review.

TIMS 2000 EN Package with Mobile Cart (no audio)

Includes 27" touch screen monitor, mobile cart, standard onsite installation & training, Video Isolator, Trigger Kit, & TIMS Review.

TIMS 2000 SP Package

Includes 24" LCD monitor, standard onsite installation & training, Video Isolator, Trigger Kit, & TIMS Review.

TIMS 2000 SP Package with Mobile Cart

Includes 27" touch screen monitor, mobile cart, standard onsite installation & training, Video Isolator, Trigger Kit, & TIMS Review.



978-458-4624
info@tims.com
www.tims.com