2653 WALK-IN BATHTUB 26" x 53"

Key Attributes

Walk-in Type:	.Walk-in	
Material:	. Acrylic	
Soaker Weight:	. 110 lbs	
Jetted Weight:	. 145 lbs	
Average Fill Capacity*:	.60-70 gal	
*Estimated with occupant inside. May vary based on occupant's volume.		

Required Electrical	Service 110v
Air Bath or Whirlpool	. 1x 15 amp GFCI
Whirlpool w/ Inline Heater	. 1x 20 amp GFCI
Whirlpool & Air Jetted	. 1x 20 amp GFCI
Soaker	.n/a


Available Therapy Systems

2653**S	Non-Jetted Soaking
2653**H	Whirlpool Jetted
2653**A	Air Jetted
2653**D	Whirlpool & Air Jetted

Key Measurements (inches)

•	-
Overall Outer:	52.875 x 25.5 x 38
Interior Width:	19
Standing Area:	22 x 15
Seat Width:	17
Seat Height:	17
Soaking Depth:	30
Drain to Side Wall:	9
Drain to Back Wall:	12
Step-in Height:	7
Door Clearance:	10.5
Faucet Deck:	5
Tile Flange:	1.2
Extension Kit:	7

It is vital that tub is completely level with all leveling feet touching the floor for the door, seal and drain to work properly.



.

- IMPORTANT! TEST THE TUB'S DRAIN AND JETS (if any) BEFORE PERMANENT
 INSTALLATION!
- **INSPECT THE TUB WHEN RECEIVED!** Report any damages or missing parts within 5 business days of delivery date.
- Before ordering or installing, please consider future access to equipment.
- Please measure tub on-site before building tub surrounding!
- Jet, pump & blower locations are predetermined and cannot be changed.
- Tub door must remain closed when not in use.

- Please note that installer(s) must be certified or warranty will be voided.
- Measurements are subject to a +/- .5" variance and subject to change without notice. Measurements here in supersede all other published specs prior to publication date shown below.
- Inline heater (if included) is designed to maintain water temperature during operation. It will not heat cold water to desired bathing temperature.
- Product must be registered for warranty activation. See user manual or visit SWCORP.com/register for more information.