

1. Part Number: 8063H-300S

2. Photograph:



3. Product Information:

| | |
|---|---|
| Description: | SteriWare® Powder Lance - 300mm – Sterile |
| Use: | Single use, disposable powder sampler |
| Material of Construction: | HDPE (High Density Polyethylene) - virgin <ul style="list-style-type: none">• Conforms to FDA CFR 177.1520• Conforms to EU Regulations 10/2011• Conforms to EC Regulations 1935/2004• BSE/TSE free |
| Method of Construction: | Ultrasonic welding |
| Nominal Length: | 401mm |
| Blade Length: | 300mm |
| Diameter of Blade: | 22mm |
| Nominal Weight of Sampler: | 63g |
| Moulding & Primary Packing Environment: | Class 100,000 Cleanroom |
| Individually Bagged? | Yes (heat sealed PE/PET bag) |
| Method of Sterilisation: | Gamma Irradiation (25-45kGy) |
| Quantity per Box: | 20 |
| Recommended Storage Conditions: | Dry and ambient temperature |
| Shelf Life: | 5 years from the date of manufacture |

4. Change Control: Sampling Systems operates a policy of continual improvement. To be kept informed of all major changes to any product in the SteriWare range please email your contact details to change.control@sampling.com

5. Country of Origin: Made and packed in the UK.

6. Repeated Use:

This is a single-use, disposable sampler. Repeated use will lead to failure.

7. Approval:

Signed:



Name: Kuldeep Flora

Position: Quality Manager

Date: 29-Jun-21