

## Specification Number: 8111H-200S

Revision: Date Issued:

02 1-Mar-21

1. Part Number:

8111H-200S

2. Photograph:



## 3. Product Information:

Description:	SteriWare Powder Funnel 200mm - Sterile
Use:	Single use, disposable powder spatula & cover
Material of Construction:	<ul> <li>High Density Polyethylene (HDPE) - virgin</li> <li>FDA acceptable CFR 177.1520</li> <li>EU Regulations 10/2011 compliant</li> <li>EC Regulations 1935/2004 compliant</li> <li>BSE/TSE free</li> </ul>
Method of Construction:	Single piece, injection moulded
Max. Diameter:	200mm
Overall Height:	165mm
OD at Bottom of Spout:	43mm
ID at Bottom of Spout:	40mm
Nominal Weight:	86g
Moulding & Primary Packing Environment:	Class 100,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	Gamma Irradiation (25-45kGy)
Quantity per Box:	10



## 4-5 Forge Mills Park, Coleshill. Warwickshire. B46 1JH. UK.

Tel: +44(0)1675 466992 Email: info@sampling.com Web: www.sampling.com

SteriWare is a registered Trademark belonging to Sampling Systems Ltd.



**Specification Sheet** 

Specification Number: Revision: Date Issued: 8111H-200S 02 1-Mar-21

Recommended Storage Conditions: Shelf Life: ATEX:

ATEX:

4. Change Control:

5 years from the date of manufacture

Dry and ambient temperature

This product has not been ATEX rated – Contact Sampling Systems for details

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to *change.control@sampling.com* 

- 5. Country of Origin:
- 6. Repeated Use:

7. Approval:

Made and packed in the UK.

This is a single-use, disposable sampler. Repeated use will lead to failure.

Signed:

Name: Kuldeep Flora Position: Quality Manager Date: 1-Mar-21



4-5 Forge Mills Park, Coleshill. Warwickshire. B46 1JH. UK. Tel: +44(0)1675 466992 Email: info@sampling.com Web: www.sampling.com