

Specification Sheet

Specification Number: Revision: Date Issued: 8890P-500S 02 21-May-19

1. Part Number:

8890P-500S

2. Photograph:



3. Product Information:

Method of Construction:

Description:	SteriWare® Powder Bottle 500ml - Sterile
Use:	Single use, disposable bottle
Material of Construction:	Polypropylene (PP) – opaque • Conforms to FDA CFR 177.1520

- Conforms to EU Regulation 10/2011Conforms to EC Regulation 1935/2004
- BSE/TSE free

Single piece, injection moulded

Nominal Volume:	500ml
Brimful Volume:	540ml
Method of Construction:	Single piece, injection moulded
Height:	106mm
Lid Diameter:	94mm
Nominal Weight:	83g
Moulding & Primary Packing Environment:	Class 100,000 Cleanroom
Individually Bagged?	Yes (heat sealed PE/PET bag)
Method of Sterilisation:	Gamma Irradiation (25-45kGy)
Number of Bottles per Box:	37
Recommended Storage Conditions:	Dry and ambient temperature
Shelf Life:	5 years from the date of manufacture



Sampling Systems Ltd.

4-5 Forge Mills Park, Coleshill. Warwickshire. B46 1JH. UK. Tel: +44(0)1675 466992 Email: info@sampling.com Web: www.sampling.com



Specification Sheet Specification Number: Revision:

Date Issued:

8890P-500S 02 21-<u>May-19</u>

ATEX:

4. Change Control:

This product has not been ATEX tested.

Sampling Systems operates a policy of continual improvement.

To be kept informed of all major changes to any product in the SteriWare range please email your contact details to change.control@sampling.com

This is a single-use, disposable device. Repeated use will lead to

- 5. Country of Origin:
- 6. Repeated Use:
- 7. Approval:

Signed:

failure.

Made and packed in the UK.

Name: Kuldeep Flora Position: **Quality Manager** 21-May-19 Date:



Sampling Systems Ltd. 4-5 Forge Mills Park, Coleshill. Warwickshire. B46 1JH. UK. Tel: +44(0)1675 466992 Email: info@sampling.com Web: www.sampling.com